


CHROMABOND® HR-X

Technical data

Hydrophobic polystyrene-divinylbenzene copolymer (PS/DVB)

SPE mode:	Reversed phase
Interactions:	Hydrophobic and π - π
Particle shape:	Spherical
pH stability:	1–14
Particle size:	85 μ m and 45 μ m
Pore size:	55–60 Å
Specific surface:	1000 m ² /g
RP capacity:	390 mg/g (caffeine in water)



Recommended application

- Pharmaceuticals / active ingredients from tablets, creams and water
- Drugs and pharmaceuticals from urine, blood, serum and plasma
- Trace analysis of pesticides, herbicides, phenols, PAH and PCBs from water

Standard protocol for CHROMABOND® HR-X

MN Appl. No. 304310



Column type:
CHROMABOND® HR-X / 3 mL / 200 mg, REF 730931

Sample pretreatment:

Individual sample preparation in reference to the compounds and matrix (adjust pH value if necessary).

Conditioning: 5 mL methanol, then 5 mL water
(do not let run the column dry!)

Sample aspiration: The prepared sample is passed through the column by vacuum or pressure (max. 1000 mL sample volume)

Washing: 5 mL water / methanol (95:5, v/v)

Drying: With nitrogen or air

Elution: 3 x 2 mL methanol

Further analysis:

Evaporation and reconstitution (if necessary); HPLC or GC

These conditions are a starting point for SPE method development. Further optimization may be required to improve results.

Good to know

A possible replacement for:

- Nexus
- ENVI-Chrom P
- Bakerbond H₂O-phobic DVB
- Strata™-X



Modern polymeric CHROMABOND® SPE phases

Applications

Determination of pyrrolizidine alkaloids


MN Appl. No. 306620

Chromatographic conditions

Columns: CHROMABOND® HR-X/85 µm/3 mL/200 mg
 MN REF: 730921
 Pretreatment: The following analysis were performed with standard solutions
 Conditioning: 5 mL methanol, 5 mL water
 Application: 10 mL neutralized standard solution with a flow rate of 3 mL/min
 Washing: 2 x 5 mL of water with a flow rate of 3 mL/min
 Drying: 5–10 min with vacuum
 Elution: 5 mL methanol
 Eluent exchange: Add 1.0 mL water as keeper. Evaporate eluate to a volume of 0.5 mL at 40 °C under a stream of nitrogen and fill up to 1.0 mL with water / methanol (95:5, v/v).

Further analysis:

HPLC determination of recovery rates with EC 150/2 NUCLEOSHELL® RP 18plus, 2.7 µm (REF 763236.20) in reference to MN Appl. No. 127480



Superior to silica based RP phase

CHROMABOND® HR-X shows higher recovery rates for most tested pyrrolizidine alkaloids than CHROMABOND® C18 ec under the given conditions.



Enrichment of opiates

MN Appl. No. 306710

Chromatographic conditions

Columns: CHROMABOND® HR-X/45 µm/3 mL/60 mg
 MN REF: 730936P45
 Pretreatment: 400 µL methanolic standard solution were diluted with 50 mmol/L phosphate buffer pH 7.0 to 20 mL. 2.5 mL of this solution are equal to 5 ng of each analyte
 Conditioning: 3 x 1 mL methanol, 3 x 1 mL water, then 3 x 1 mL 50 mmol/L phosphate buffer pH 7.0
 Aspiration: 2.5 mL of pretreated sample solution is passed through the column at a flow of 1–2 mL/min
 Washing: 3 x 1 mL 50 mmol/L phosphate buffer pH 7.0, 3 x 1 mL water
 Drying: 5 mL air by pushing with a syringe
 Elution: 3 x 1 mL 0.1 % formic acid in methanol

Solvent change: Eluate is evaporated to dryness at 30 °C under a stream of nitrogen and then redissolved in organic solvent suited for the subsequent analysis.

Further analysis:

HPLC determination of recovery rates with EC 100/2 NUCLEOSHELL® Biphenyl, 2.7 µm (REF 763634.20) in reference to MN Appl. No. 128880

Compound	Recovery rate [%]	Standard deviation [%]
Ecgonine methyl ester	94	0
Morphine	77	3
Dihydrocodeine	101	1
Codeine	97	1
6-Acetylmorphine	89	1
Benzoylcegonine	102	0
6-Acetylcodeine	100	0
Cocaine	109	1
Noscapine	95	1
Papaverine	98	2

