


Operating Manual

Micropump mp6 and Controller



Last update: 2010-08-25



Table of contents

1. General 4

 1.1 Declaration of conformity4

 1.2 Description of functions4

2. Proper use 5

 2.1 Intended purpose5

 2.2 Misuse5

 2.3 Staff selection and qualification.....5

 2.4 About this operating manual5

3. Technical specifications 6

 3.1 Technical specifications mp6 ^{1,2}6

 3.2 Typical flow characteristics7

 3.3 Flow curves for selected media8

 3.3.1 Medium: DI water.....8

 3.3.2 Medium: air.....9

 3.4 Final inspection9

4. Operating the micropump..... 9

 4.1 Connecting the micropump.....9

 4.1.1 Connecting the mp6 to the cable 10

 4.2 Cleaning the system..... 11

5. The mp-x controller 12

 5.1 Technical specifications mp-x controller 12

 5.2 Electrical signal form 13

 5.3 Connecting the pump to the mp-x 14

 5.4 Operation of the mp-x 15

 5.4.1 Installation of the drivers "USB Micropump Control" and "USB Serial Port" 16

 5.4.2 Operation via USB port (after installation of the drivers)..... 16





- 6. The mp6-a controller 17**
 - 6.1 Technical specifications mp6-a 17
 - 6.2 Connecting the mp6 to the mp6-a 18
 - 6.3 Operation of the mp6-a with batteries 18
 - 6.4 Operation of the mp6-a with mains adaptor 18
- 7. The mp6-OEM controller 19**
 - 7.1 Technical specifications mp6-OEM 19
 - 7.2 Typical flow characteristics 19
 - 7.3 Connecting the mp6 to the mp6-OEM controller 20
 - 7.4 Typical schematics 20
 - 7.4.1 Fixed pump rate 20
 - 7.4.2 Pump rate set by external components 20
 - 7.4.3 Pump rate control via microcontroller 22
 - 7.5 Electrical characteristics 22
 - 7.6 Pin description 23
 - 7.7 Noise reduction 23
- 8. The mp6-OEM-EVA evaluation board 24**
 - 8.1 Safety notice 24
 - 8.2 Electrical specifications mp6-OEM-EVA evaluation boards 24
 - 8.3 Functional elements 24
 - 8.4 Operation 25
 - 8.4.1 Pump frequency setting with jumper J1 25
 - 8.4.2 Pump amplitude setting with jumper J2 25
 - 8.4.3 Operation voltage setting with jumper J3 25
 - 8.5 Connecting the pump via CON2 26
- 9. Passive check valve mp-cv 27**
- 10. Trouble shooting 28**



1. General

This operating manual contains all necessary instructions for the installation, commissioning, operation and maintenance of the mp6 micropump and of the controller mp-x, mp6-a and mp6-OEM with the mp6.

The unit has been designed with state-of-the-art technology and in accordance with all relevant safety regulations. However, a risk of damage to the units, other property, the operator and/or other persons cannot be fully excluded.

Always comply with the following general instructions:

- Before working with a pump, you must be fully familiar with its operation and functions.
- Prior to operating the pump, read this operating manual and adhere to all instructions.
- Refrain from any operations that might endanger the safety of the unit.

Bartels Mikrotechnik GmbH rejects any responsibility for damages to persons or property resulting from non-compliance with the instructions in this manual. In this case all warranties shall be void.

1.1 Declaration of conformity

Bartels Mikrotechnik GmbH is certified according to DIN EN ISO 9001:2000 and declares that the products are compliant to the RoHS directive and the controller comply with the requirements of EMC 89/336/EEC and CE markings have been affixed to the devices.

1.2 Description of functions

The micropumps have been developed for the transport of gases or liquids. The mp-x controller, the mp6-a controller and the mp6-OEM controller have been developed for operating one mp6.

Bartels Mikrotechnik can assume no liability for damages resulting from the pump media. This applies especially for hazardous fluids.

The pumps must be operated with Bartels Mikrotechnik electronics. Bartels Mikrotechnik GmbH cannot guarantee the proper work of the units with customer specific electronics. If other controllers than the ones from Bartels Mikrotechnik are used, Bartels Mikrotechnik disclaims any warranty.

Please make sure that only skilled personnel works with the pump control and micropump. The micropump shall be under constant supervision under running conditions. And please note that components of the controller and pump are operating with high-voltage. Therefore persons wearing pacemakers are recommended to avoid the operating system.

Bartels Mikrotechnik assumes no liability for abnormal handling, improper or negligent use of the micropump and the controller that is not conform to the specified purpose of the system. This applies especially for micropump controllers, components and systems of other manufacturers, which have not been certified by Bartels Mikrotechnik.

We guarantee that the micropumps comply to the actual state of scientific and technical knowledge and due to this the operational risks are limited to a minimum.

Do not open the housing of the micropump and the controllers. In those cases Bartels Mikrotechnik can not issue a guaranty anymore. Please keep this manual safe and give a copy to all users.

2. Proper use

2.1 Intended purpose

The micropump is intended for pumping liquids or gases with varying flow rates controlled by the electronics. The mp-x, the mp6-a and the mp6-OEM controllers are intended for operating one mp6. Any other use of the micropump or controller unit is deemed improper.

Do not make any modifications or extensions to the pump or controller without the prior written consent of the manufacturer. Such modifications may impair the safety of the unit and are prohibited. Bartels Mikrotechnik GmbH rejects any responsibility for damage to the unit caused by unauthorized modifications to the pump and risk and liability are automatically transferred to the operator.

2.2 Misuse

The use of liquids, which may alone or in combination create explosive or otherwise health-endangering conditions (including vapors) is not permitted.

2.3 Staff selection and qualification

All work in connection with the installation, assembly, commissioning/decommissioning, disassembly, operation, servicing, cleaning and repairing of the pump and the controller must be carried out by qualified, suitably trained and instructed personnel.

Work on electrical components and assemblies must be carried out by personnel with the necessary qualifications and skills.

2.4 About this operating manual

Warnings and important notes are clearly identified as such in the text. The relevant text sections feature a specific sign. However, this icon cannot replace the safety instructions. Therefore, carefully read all safety instructions in this manual. Warnings and important notes in this text are highlighted as shown below, according to the severity of the damage that might result from non-compliance.



⚠ DANGER
DANGER INDICATES A HAZARD WITH A HIGH LEVEL OF RISK WHICH, IF NOT AVOIDED, WILL RESULT IN DEATH OR SERIOUS INJURY.



3. Technical specifications

3.1 Technical specifications mp6 ^{1,2}

Pump type	piezoelectric diaphragm pump		
Number of actuators	2		
Dimensions without connectors	30 x 15 x 3.8 mm ³		
Weight	2 g		
Fluidic connectors	tube clip (outer diameter 1,6 mm, length 3,5 mm)		
Electric connector	flex connector / Molex FCC 1,25 mm pitch		
Power consumption	< 200 mW		
Self-priming	yes ³		
Pumped media	liquids, gases and mixtures		
Operating temperature	0 – 70°C ⁴		
Life time	5000 h ⁴		
IP code	IP33 ⁵		
Materials in contact with media	polyphenylene sulphone (PPSU)		
Suitable evaluation controller	mp-x and mp6-a		
Typical values of flow and back pressure for selected media (values defined with mp-x: 250 V, SRS):			
Gases	Max. flow	on request	
	Max. back pressure	on request	
Liquids	Water	Max. flow	6 ml/min +/- 15% (100 Hz)
		Max. back pressure	550 mbar +/- 15% (100 Hz)

¹ Preliminary data

² Content is subject to changes without notice.


³ Conditions: suction pressure < 10 mbar, DI water, settings mp-x: 100 Hz, 250 V, SRS, the max. flow rate will be reached after a few minutes of operation time.

⁴ Value of previous version

⁵ Can be changed to IP44



3.2 Typical flow characteristics





The flow rate of the pumps shows a linear dependency on the back pressure.

At 0 mbar back pressure the maximum pump rate can be achieved and at the maximum back pressure the flow rate is decreased to 0 ml/min.

Increase in amplitude linearly increases the flow rate to the maximum.


The amplitude defines the stroke of the actuator and therefore the displacement of the pumping media per pump cycles.



Rising frequencies result in a maximum flow rate at resonance frequency. At frequencies above the resonant point the flow rate decreases again.



The resonance frequency and the maximum flow rate strongly depend on the viscosity of the media.

The lower the viscosity, the higher the maximum flow rate and the resonance frequency.



3.3 Flow curves for selected media

3.3.1 Medium: DI water



3.3.2 Medium: air

Data is not available yet, please contact us for more information.


3.4 Final inspection

After production the micropumps have to pass a final inspection. They are tested concerning the maximum flow of 6 ml/min and maximum back pressure of 550 mbar (DI water; settings mp-x: 100 Hz, 250 V, SRS) and the self-priming characteristic (conditions: suction pressure < 10 mbar; DI water; settings mp-x: 100 Hz, 250 V, SRS, the max. flow rate will be reached after a few minutes of operation time).

In order to guarantee proper function of the delivered goods and exclude transportation damages please check the incoming devices according to specifications after receipt. On the basis of these results a replacement can be carried out within 14 days after delivery free of costs.

4. Operating the micropump

4.1 Connecting the micropump



Please connect suitable tubes to the inlet and outlet. The tubing should have an inner diameter of 1.3 mm. The micropump has to be connect to a suitable controller, described in the next sections.

In general the micropump can be driven with positive alternating voltages with maximum amplitude of 250 V at a frequency between 0 and 300 Hz. A rectangular signal results in best fluidic performance while a sine wave minimizes the audible noise. The actuators must be driven with a 180° phase shift in the signal.


If a pump gets damaged while using a customer's controller we do not provide any warranty.

! DANGER

THE MP6-SERIES MICROPUMPS ARE OPERATED AT HIGH VOLTAGES. PEOPLE WHO WEAR HEART PACE-MAKERS ARE THEREFORE STRONGLY ADVISED NOT TO COME NEAR THEM.

4.1.1 Connecting the mp6 to the cable


The mp6 pump can be connected via a standard electrical connector manufactured by Molex. The layout of the connector and pin assignment of the mp6 is as shown below. Each piezo (P1 / P2) has a single lead for the negative (-P1/-P2) and the positive (+P1/+P2) supply voltage.



The recommended connector is a 4 pin 1.25 mm pitch FCC connector manufactured by Molex, part number 39 53 2045. The connector is available in different versions to fit various applications. Reference drawings and further specifications are available under www.molex.com.

The recommended maximum wire length between controller and the pump is 1 m. The customer has to assure that the assembly of pump, controller and electrical connection complies with the EMC regulations and electrical safety in the specific field of application.

To connect the mp6 and Molex, refer to following three figures in the right column. Orientate both components as indicated on the first picture, the mp6 with the metallic contacts down. Then insert the mp6 flex into the Molex connector. Close the Molex connector to complete the interconnection between both components.



In order to prevent damage to the flexible cable, the following points should be considered for the final pump assembly

- the flexible connector must not be bent around sharp edges or kinked
- the flexible connector must not be bend on top or bottom of the pump
- the Molex connector is not water tight, additional sealing with e.g. silicone might be necessary
- fixation of the Molex connector in the final assembly is recommended



4.2 Cleaning the system

The pump can be washed with water, alcohol (iso-propanol) or if necessary with weak acid by pumping or by flushing with the help of a syringe.

!! Only use the syringe pulling from the outlet in pumping direction!!



5. The mp-x controller


5.1 Technical specifications mp-x controller

Dimensions	7.5 x 16 x 20 cm ³
Weight	approx. 800 g
Pumping media	liquids or gases
Max. flow	6 ml/min (water)
Controlling parameters	amplitude, frequency, signal form
Amplitude range	0 – 250 V
Frequency range	0 – 300 Hz
Signal form	SRS, rectangular, sine
Power supply	mains adaptor
Current consumption	500 mA with 7.5 V
USB (port)	one, a CD with driver software is enclosed.



5.2 Electrical signal form

SRS:




mpx, SRS-Signal, 250Vpp, 100Hz

Rectangular:



mpx, Rectangel-Signal, 250Vpp, 100Hz


Sine:



mpx, Sine-Signal, 250Vpp, 100Hz



5.3 Connecting the pump to the mp-x



- 1 Main power switch
- 2 Power supply connector
- 3 USB interface for computer connection
- 4 Two connectors for the mp5
- 5 One connector for the mp6

Please note that it is only possible to connect either one mp6 or a maximum of two mp5 to the mp-x otherwise a maximum voltage drop is possible!

Step 1: Plug the micropump control cable into the corresponding micropump connector.

Step 2: Check the mains adaptor plug polarity. It is pictured next to the power supply connector at the back of the controller. If the plug polarity is wrong, the controller can't work. Please make sure that the setting on the included connector is attuned to 7.5 V.

Step 4: Connect the mains adaptor with the power supply connection.

Step 5: Plug the mains adaptor into a mains socket.

Step 6: Now you can start the control unit with the main power switch.

⚠ DANGER


**THE "MICROPUMP OUT" CONNECTOR CAN CARRY HIGH VOLTAGE !
BE CAREFUL, WHEN YOU PLUG IN THE MICROPUMP CONTROL CABLE!**

5.4 Operation of the mp-x

The mp-x provides three parameters to be selected independently of each other to control the micropumps:

Frequency, amplitude and signal wave form.

It is possible to change the settings while the pump is operating. But to extend the life time it is advisable to turn off the micropump first before changing the settings.



- 1 LCD display
- 2 On/Off switch: push the on/off switch to turn the pump and control diode on or off
- 3 Save: push the save switch to save the present settings
- 4 Frequency (0 – 300 Hz): push the switch up to raise and down to reduce the frequency
- 5 Signal form modus: push the switch once to choose between SRS, rectangular or sine
- 6 Amplitude (0 – 250 V): push the switch up to raise and down to reduce the amplitude

To operate the micropump, prepare the controller as described in point 5.3 and follow the steps below:

Step 1: Choose a frequency by pushing the frequency switch up and down.

Step 2: Choose one of the signal forms by pushing the signal switch.

Step 3: Set the amplitude by pushing the amplitude switch up and down.

Step 4: Push the on/off switch to turn on the micropump and the control diode.

Step 5: Push the on/off switch again and the micropump will stop while the control diode turns off.

Before the control unit is turned off, you can save the present settings by pushing the save switch.

For shutdown of the pump please switch off the mp-x first and then disconnect the power plug. Do not unplug the micropump before switching off the controller.

5.4.1 Installation of the drivers "USB Micropump Control" and "USB Serial Port"

Step 1: Connect the control unit with the USB port and turn it on. A message appears that new hardware was found and the hardware assistant starts automatically. Please click "next" to continue.

Step 2: Choose "find a suitable driver for the device" and click "next" to continue.

Step 3: Activate the box "search for new driver" and insert the accompanying CD.

Step 4: After you have found the driver "USB Micropump Control", click "next" to start the installation. Maybe a message appears that the windows-logo-test has failed. This warning can be ignored. After installation click "complete" to finish the installation.

Step 5: After finishing the installation the assistant will start again to install the corresponding USB serial port by repeating the described above.

5.4.2 Operation via USB port (after installation of the drivers)

The driving parameters can be set via your PC. For this purpose you can use any software (or programming language) capable of sending commands to a COM-Port, like the Hyperterminal shown in the example below. As the pump will be switched off during data transfer, this is not suitable for continuous control.

Example with Hyperterminal

Step 1: Connect the control unit to your computer and turn it on.

Step 2: Start Windows Hyperterminal. Every new session has to be titled.

Step 3: Choose the com-port specified in the device manager.

Step 4: The connection-settings have to be (9600, 8, n,1).

Possible commands (followed by the enter key)

bon	turns the micropump on
boff	turns the micropump off
<i>F(1-300)</i>	sets the required frequency between 1 and 300 Hertz:
F100	...for example 100 Hertz.
<i>A(0-250)</i>	sets the required amplitude between 5 and 250 Volt:
A100,5	...for example 100.5 Volt. (in increments of 0.5 Volt)
MS	sets signal form modus (S)inus
MR	sets signal form modus (R)ectangle



MC	sets signal form modus SRS
(enter key)	displays present settings of the control unit

If required we can send you a package of LabView-Routines as an example of implementation.

6. The mp6-a controller


6.1 Technical specifications mp6-a

Dimensions	7 x 5 x 2 cm ³
Weight	approx. 80 g (including batteries)
Pumping media	liquids
Max. flow	4.5 ml/min (water)
Controlling parameters	frequency
Frequency range	25 - 120 Hz
Amplitude	230 Vpp
Signal form	similar to rectangular
Power supply	2 AA batteries ¹ or mains adaptor
Current consumption	approx. 50 mA (3 V)



¹ The flow range decreases by potential drop of battery.

6.2 Connecting the mp6 to the mp6-a



Step 1: Connect the mp6 with the control unit. For connecting the pump to the cable see chapter 4.1.1.

Step 2: Adjust the flow rate by turning the potentiometer.

6.3 Operation of the mp6-a with batteries

Slide the lid off, insert two 1.5 V AA batteries and switch on the micropump with the On/Off – switch. The flow range decreases by potential drop of battery.

6.4 Operation of the mp6-a with mains adaptor

Connect the controller with the included adapter, plug it in a mains socket and switch on the micropump with the On/Off – switch and it will operate immediately. Please make sure that the setting on the included connector is attuned to 3 V. The function of the controller and the micropumps will not be influenced, if the batteries have not been removed.




7. The mp6-OEM controller

The mp6-OEM is a small, easy to use low cost driving circuit developed for the mp6 micropump. It generates up to 235 V peak to peak voltages from a 3-5 V supply.

Its low power consumption makes it ideal for battery powered handheld devices or even solar powered devices. The module can be integrated into a PCB design like a 14 pin DIL package; due to its resin encapsulation it is robust in usage.


Build in interface allows the user to adapt frequency and/or amplitude to its application by the use of a few additional components or a microcontroller.



7.1 Technical specifications mp6-OEM


Dimensions	12 x 22 x 10 mm ³
Pumping media	liquids or gases
Adjustable parameters	amplitude / frequency
Amplitude range	100 – 235 V
Frequency range	25 – 120 Hz
Signal form	Similar to rectangular
Power supply	2.5 – 5.5 V DC
Current consumption	Approx. 40 mA at 3 V
Max. flow rate mp6 (typ.)	4.5 ml/min (water)

7.2 Typical flow characteristics



7.3 Connecting the mp6 to the mp6-OEM controller


The mp6 pump can be connected to the mp6-OEM via a standard electrical connector manufactured by Molex, see chapter 4.1.1

 DANGER
THE MP6-OEM MODULE CAN CARRY HIGH VOLTAGE ! BE CAREFUL, WHILE CONNECTING AND HANDLING THE MODULE!

7.4 Typical schematics

7.4.1 Fixed pump rate

The mp6-OEM can drive the mp6 without the need of external components. In this case the pumps frequency and amplitude is fixed determined by internal components to 235 V and 100 Hz typical.




Schematic 1: Fixed amplitude of 235 V and fixed frequency of 100 Hz $\pm 10\%$ with internal components.

7.4.2 Pump rate set by external components

In this example the frequency and the amplitude is set by external components. The amplitude can be varied from 100 V to 235 V peak to peak. The frequency can be changed from 25 Hz to 120 Hz. For higher frequencies the amplitude will decrease. Lower frequencies are also possible but stability of the output signal needs to be checked.

To set the amplitude, either a potentiometer with 10 kOhms, or a voltage divider of two resistors R1 and R2 as shown in the schematic 2 can be used. The following table shows typical resistor values. Voltages at the

AMPLITUDE pin below 0,5 V will not decrease the pump amplitude further, voltages larger than 1,3 V will not increase the pump amplitude further.



Schematic 2: Frequency and amplitude set with external components

Amplitude of micropump [V]	VDD [V]	R1 [Ohm]	R2 [Ohm]	
100V (minimum)	3V	1,8k	10k	Or amplitude pin connected to GND
200 V	3V	12k	22k	
~ 235 V (maximum)	3V	12k	33k	Or amplitude pin connected to VDD
100V (minimum)	5V	1k	10k	Or amplitude pin connected to GND
200 V	5V	8,2k	33k	
~ 235 V (maximum)	5V	12k	33k	Or amplitude pin connected to VDD

Table : Typical resistor values for different pump amplitude values

To set the frequency, a capacitor C1 between 1 and 10 nF can be used as shown in the schematic above. Typical capacitor values are shown in the following diagram.





Table : Typical capacitor values for different pump frequency values

7.4.3 Pump rate control via microcontroller

To control the mp6-OEM with a microcontroller, an external clock with four times the desired frequency connected to the CLOCK pin can directly set the output frequency. The amplitude can either be adjusted by an analog voltage between 0.5 V and 1.3 V (higher voltages will not increase the amplitude), or by an equivalent PWM signal with a frequency between 0.2 and 1 MHz connected to the AMPLITUDE pin. The PWM signal should be higher than 1.3 V to make the output solely dependent on the duty cycle of the PWM signal.

For a minimum of power consumption the electronics can be switched off by applying zero volts to the AMPLITUDE input and the SHUTDOWN port, but the most efficient way is to cut the power supply.



Schematic 3: External control via microcontroller

7.5 Electrical characteristics

One mp6 connected, internally defined frequency and slew rate

Parameter	Symbol	Conditions	Min	Typ	Max	Units
Power Supply Voltage	VDD		2.5		5.5	V
Power Supply Current, average	IDD	VDD = 5 V (1)			30	mA
Control voltage AMPLITUDE			0,5		1.3	V
Peak to Peak Output Voltage	V _{pump}	AMPLITUDE = 1.3 V	230	235	245	V
Peak to Peak Output Voltage	V _{pump}	AMPLITUDE = 0.5 V	90		100	V
PWM frequency AMPLITUDE			0.2		1	MHz
Internal Pump Frequency	F	VDD = 5 V (1)	90	100	110	Hz
Digital Inputs Low				0		V
Digital Inputs High			2V		VDD	V
Capacity at Clock Input			1.0	2.2	10	nF
Input Current AMPLITUDE			1		3	μA
Current in shutdown mode				1,6		μA

(1) Output signal set by internal components

7.6 Pin description

VDD	Power supply voltage
GND	Ground
SHUTDOWN	To shutdown the device, AMPLITUDE <u>and</u> SHUTDOWN needs to be tied to GND.
CLOCK	Output frequency control, the frequency can be set to nominal 100Hz by connecting this pin to CLOCK_INT (Schematic 1) A capacitor of 1 to 10 nF can be connected between this pin and GND to set another frequency (Schematic 2) The output frequency can be set by a clock signal with <u>four times</u> the desired output frequency
CLOCK_INT	Output frequency control, the frequency can be set to nominal 100Hz by connecting this pin to CLOCK
AMPLITUDE	Apply a DC Voltage (0-1.3 V) or a PWM signal (0.2-1 MHz) to this input to adjust the amplitude of the output from 100 V to 235 V
SLEWRATE SLEW_INT	Slew rate control. This pin is connected to SLEW_INT Internal slew rate resistance, connect this pin to SLEW
GND2, GND3	Internally connected to GND, can be left unconnected
+P1	Piezo 1 positive (see connection diagram for the mp6)
-P1	Piezo 1 negative (see connection diagram for the mp6)
+P2	Piezo 2 positive (see connection diagram for the mp6)
-P2	Piezo 2 negative (see connection diagram for the mp6)

7.7 Noise reduction

If the noise generated by the pump is critical a series resistor of 2-10 kΩ in the P1+ and the P2+ line between the MP6-OEM and the pump will help. There is no limit for the resistor value but it will decrease the maximum pump performance.



8. The mp6-OEM-EVA evaluation board

8.1 Safety notice


The mp6-OEM generates voltages of up to 250 V peak to peak. All parts of the mp6-OEM evaluation board can carry voltages in this range. Therefore the board should only be used by qualified personal. Although the output power of the module is very low, proper insulation according to the application conditions needs to be considered by the customer. This especially applies to the lower side of the PCB. Contact with water or other liquids needs to be prevented. The pump must not be changed while a driving voltage is applied to the board.



8.2 Electrical specifications mp6-OEM-EVA evaluation boards

As the evaluation board is based on the mp6-OEM module, all electrical characteristics and specifications of this product must be considered. Please see chapter 7 of this manual for more detail.

8.3 Functional elements



Elements are listed with their names according to the printed description on the PCB

Connectors:

- CON 1 – Screw terminal for external power supply and external clock / amplitude signal
- CON 2 – Solder terminal for extension cable to connect one mp6 micropump
- CON 3 – Molex connector to connect one mp6 micropump
- USB connector for voltage supply via USB

Jumpers:

- JP1 – Jumper for pump frequency setting
- JP2 – Jumper for pump amplitude setting
- JP3 – Jumper for power supply setting

Others:

- P1 – Variable resistor for amplitude adjustment



8.4 Operation

To operate a pump with the evaluation board, the following steps are necessary:

- 1) Connect the mp6 pump to the board according to the description in chapter 4.1.1. Due to the orientation of the connector, the pump needs to be inserted with its metallic contacts upwards.
- 2) Choose the pump frequency setting with Jumper 1
- 3) Choose the pump amplitude setting with Jumper 2
- 4) Choose the power supply setting with Jumper 3
- 5) Connect the board with the voltage source

8.4.1 Pump frequency setting with jumper J1


Setting of Jumper J1



Internal frequency of the mp6-OEM (100 Hz)



Frequency defined by capacitor C1 (predefined to 50 Hz)




Frequency defined by CLK input on terminal CON1 – Pin 3

More information on the frequency setting with the CLK signal can be found in chapter 7.4.2. The capacitor C1 can as well be changed according to chapter 7.4.2 by resoldering.

8.4.2 Pump amplitude setting with jumper J2


Setting of Jumper J2



Maximum voltage (235 V)



Amplitude defined by variable resistor P1



Amplitude defined by AMP input on terminal CON1 – Pin 4

More information on the amplitude setting with the AMP signal can be found in chapter 7.4.1.

8.4.3 Operation voltage setting with jumper J3

Setting of Jumper J3



Driving voltage via screw terminal CON 1
Pin 1 (Vcc) and Pin2 (GND)



Driving voltage via USB port

8.5 Connecting the pump via CON2

If the mp6 micropump should not be connected directly to the PCB with the Molex connector CON 3, an extension cable can be soldered to the connector CON 2.

The solder pads have the following pin assignment:

■	P2 +	(positive voltage piezo 2)
■	P2 -	(negative voltage piezo 2)
■	P1 -	(negative voltage piezo 1)
■	P1 +	(positive voltage piezo 1)

CON 2

For further details, please refer to chapter 4.1.1 of this manual. Please make sure that the cable can handle voltages up to 250 V peak-to-peak and ensure proper insulation of the cable.




9. Passive check valve mp-cv

When the micropump is switched off, a back flow of the pumping medium, depending on a differential pressure between in- and outlet is possible. In order to impede a back flow Bartels Mikrotechnik offers a passive check valve integrated in stainless steel.

The connection of the valves with the micropump can be done easily with suitable tubing. The valve should be placed between the micropump and the outlet reservoir.

Dimensions:



The valves are individually marked regarding the free flow direction.

Technical specifications:

Material	silicone, stainless steel
Fluidic connectors	barb inner diameter: 1.6 mm length: 5.6 mm
Cracking pressure	typically < 35 mbar
Max. back pressure	500 mbar
Leak rate	< 20 μ l/h for liquids (500 mbar)

Please note that the volume flow of the micropump will be influenced by the check valve.



10. Trouble shooting

Observation	Trouble shooting
Maximum flow rate cannot be achieved	<ul style="list-style-type: none"> - Gas bubbles within the system, the compressible gas volume in the pump can result in decrease of pump rate = system should be primed with a syringe. Manual priming should be carried out by flushing the pump with the help of syringe. Only pull the fluid from the outlet in pumping direction. Do not exceed 1 bar of internal pressure. - Tubing is too long = shorten tubing - Tubing diameter is too small = use appropriate tubing with an inner diameter of 1.3 mm - Back pressure is too high = reduce pressure or lower outlet reservoir - Check signal form used to drive pump = highest flow rate is achieved by SRS signal and at resonance frequency (100 Hz for water) - The flow rate is viscosity dependent as described in 3.2
Non-linear flow behavior	<ul style="list-style-type: none"> - The linear range of the pump characteristic is described for selected media in 3.2 - Gas bubbles in the pump act as compressible volume and when they are digested through the pump this can lead to non-linear flow behavior = system should be primed with a syringe as mentioned above
Fluid is flowing through the system although pump is switched off	<ul style="list-style-type: none"> - The valves inside the pump are opening and closing by differential pressure of in- and outlet = lower in- or outlet reservoir to avoid hydrodynamic behavior or add check valve as described in 8.
No buzzing sound = pump is not working	<ul style="list-style-type: none"> - Please check if the mains adaptor is correctly connected - Defect piezo actuator = contact Bartels Mikrotechnik for support





Bartels ikrotechnik

Bartels Mikrotechnik GmbH

Emil-Figge-Str. 76a, 44227 Dortmund
Germany

www.bartels-mikrotechnik.de

microComponents@bartels-mikrotechnik.de

Tel: +49-231-9742-500

Fax: +49-231-9742-501

