Use and maintenance

**REFRIGERATORS / FREEZERS**

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**KÜHLSCHRÄNKE / FROSTER**

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**REFRIGERADORES / CONGELADORES**

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**REFRIGERADORES / CONGELADORES**

**fiocchetti**  
THE COLD MANUFACTURER

**ENGLISH** **Pag. 4**

**FRANÇAIS** **Pag. 25**

**DEUTSCH** **Pag. 46**

**ESPAÑOL** **Pag. 67**

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<b>1</b>	<b>GENERAL INFORMATIONS</b>	<b>Pag.5</b>
1.1	CERTIFICATION	
1.2	TESTING AND GUARANTEE	
1.3	SCOPE,CONTENT AND ADDRESSEES OF THE MANUAL	
1.4	CLIENT'S RESPONSIBILITIES	
1.5	INSTRUCTIONS FOR SERVICE REQUEST	
<b>2</b>	<b>PRODUCT DESCRIPTION</b>	<b>Pag.6</b>
2.1	TECHNICAL DESCRIPTION	
<b>3</b>	<b>SAFETY</b>	<b>Pag.7</b>
3.1	GENERAL SAFETY REGULATIONS	
3.2	APPLICATIONS AND INTENDED USE	
3.3	USE CONTRAINDICATION	
3.4	SAFETY AND ACCIDENT PREVENTION	
3.5	SAFETY PROTECTIONS ADOPTED	
<b>4</b>	<b>TRANSPORT AND HANDLING</b>	<b>Pag.8</b>
4.1	TRANSPORT AND HANDLING	
4.2	POSITIONING	
4.3	WIRING AND ELECTRICAL HOOK-UP	
4.4	SET-UP OPERATIONS	
4.4.1	PRELIMINARY CHECKS	
4.4.2	INDICATION FOR OPTIMAL DUTY	
4.4.3	SPECIAL WARNINGS FOR MODELS ARTIC-SUPERARTIC	
4.4.4	INTERNAL FITTING NEW MODELS	
4.4.5	INTERNAL FITTING	
4.4.6	INTERNAL FITTING WITH BRACKETS	
<b>5</b>	<b>FUNCTIONING</b>	<b>Pag.13</b>
5.1	COMPLETE CONTROL PANEL	
5.2	THERMOSTAT	
5.3	VISUAL AND ACOUSTIC ALARM SYSTEM	
5.4	FUNCTIONING AND MAINTENANCE OF THE TEMPERATURE CHART RECORDER	
5.4.1	CHANGING OF THE PAPER CHART	
5.4.2	REPLACING OF THE BATTERY	
5.4.3	REPLACING OF THE INK-TIP	
<b>6</b>	<b>ORDINARY CARE AND MAINTENANCE</b>	<b>Pag.18</b>
6.1	PROHIBITION OF THE REMOVAL OF THE SAFETY DEVICES	
6.2	CLEANING OF THE STRUCTURE, INTERNAL AND EXTERNAL	
6.3	CONDENSER CLEANING	
6.4	CONDENSATE WATER DRAINING	
<b>7</b>	<b>EXTRAORDINARY MAINTENANCE AND REPARATIONS</b>	<b>Pag.20</b>
7.1	PROTECTIONS REMOVAL	
7.1.1	BOTTOM FRONT PANEL	
7.1.2	TOP FRONT PANEL AND ELECTRIC CIRCUIT COVER	
<b>8</b>	<b>DEMOLITION</b>	<b>Pag.22</b>
<b>9</b>	<b>ATTACHMENTS</b>	<b>Pag.23</b>
<b>10</b>	<b>DATA PLATES</b>	<b>Pag.23</b>
10.1	DATA PLATES WITH THE CHARACTERISTICS OF THE APPLIANCE	
10.2	OTHER INDICATION LABELS	
<b>11</b>	<b>DIAGNOSTIC</b>	<b>Pag.24</b>

## 1.1 CERTIFICATION

All the refrigerated cabinets are manufactured in conformity to the EC directives applicable at the moment of the emission of the product on the market.

All the refrigerated cabinets are certified according to the directive CEE 73/23 and further amendments and manufactured according to the standards for the electrical appliances for laboratory application (CEI EN 61010-1).

## 1.2 TESTING AND GUARANTEE

The appliance is tested in our workshop in compliance with established regulations and then shipped ready for use.

The guarantee is valid for a full 12 months from the date of delivery of the appliance and it covers the repair or replacement of defective parts, with the exception of electrical and electronic components.

Manifest defects or difference with respect to the client's order must be communicated to the manufacturer within five days from the receipt of the goods or they will not be covered by the guarantee terms.

Any hidden or other defects must be communicated to the manufacturer within five days from the time that they are discovered and, in any event, within the maximum guarantee term of six months. The purchaser shall be entitled only to request repair or replacement of the goods. The purchaser is not entitled to claim compensation for direct or indirect damages of any whatsoever nature. In any event, the entitlement to repair or replacement of the materials must be exercised within the maximum term of the guarantee, which is contractually, stipulated to cover a shorter period than the maximum term established by law.

Repairs or replacement of defective materials will be carried out at the manufacturer's workshop; material returned to the manufacturer must be shipped carriage paid and will be returned to the purchaser carriage forward.

## 1.3 SCOPE, CONTENT AND ADDRESSEES OF THE MANUAL

This manual has been prepared with the scope of supplying all the instructions required for the correct use of the appliance and to maintain it in optimal condition. It also contains important user safety information. The following professional roles are explained in order to define the responsibilities of each:

**Installer:** qualified technician who positions the appliance and places it in service it in accordance with the instructions in this manual.

**User:** the person who, after having read this manual carefully, operates the appliance in accordance with the intended use specified in this manual. The user is obliged to read the manual attentively and refer to the information in the manual at all times.

**Routine maintenance technician:** qualified technician able to perform routine maintenance of the appliance by following the instructions in this manual.

**Special maintenance technician:** qualified technician, authorized by the manufacturer to perform extraordinary maintenance of the appliance.

The manufacturer declines any whatsoever responsibility in the case of improper use of the appliance deviating from the reasonably construed intended use, and for all operations carried out that are not in compliance with the instructions laid down in the manual.

This manual must be conserved in a place that is accessible and known to all the operators (installer, user, routine maintenance technician, special maintenance technician).

This manual must not reproduced or divulged, in whole or in part, using any whatsoever means or in any whatsoever form.

#### **1.4 CLIENT'S RESPONSIBILITIES**

The customer is required to:

- execute the electrical connection of the appliance
- prepare the place of installation
- provide consumable materials for cleaning
- perform routine maintenance

In the case of power failures or malfunctions do not open the door in order to maintain uniform temperature inside the unit. If the problem persists for more than a few hours, move the material contents to a more suitable place.

#### **1.5 INSTRUCTIONS FOR SERVICE REQUEST**

For all technical problems and any requests for technical service, refer exclusively to your local dealer (See space in the last page) or directly to the manufacturer, specifying model and serial number.

## **2**

## **PRODUCT DESCRIPTION**

### **2.1 TECHNICAL DESCRIPTION**

The appliances are of these series are of the Refrigerators and Freezer in which the production of the cold to the inside is done for vaporization at low pressure in a thermal exchanger (evaporator) of a refrigerant liquid, type HCFC or HFC; the vapour therefore obtained is bought to the original state of liquid through a mechanical compression at higher pressure (compressor), followed from a cooling stage in an other thermal exchanger (condenser).

The correct and uniform distribution of the cold air into the chamber is granted through one or more electro-mechanical fans (depending from the model).

The appliance comprises a modular single body with various materials and insulation in expanded polyurethane foam, density 43 Kg/m<sup>3</sup>.

The appliance instruments are located on the front panel that in some models closes the front of the motor until, inside which the condenser unit and electrical wiring can be housed.

The appliance interior is fitted with suitable supports for wire shelves (grids) or extractible drawers. The doors are fitted with an automatic device and magnetic seal elements easily replaceable.

During the design and construction stages, all the measure have been adopted to implement total safety including radiused interior corners, funnel-shaped base panel to convey condensate to exterior, no rough surface, fixed guards protecting moving or potentially dangerous parts.

### 3.1 GENERAL SAFETY REGULATIONS

Read this manual carefully and follow the prescriptions contained herein.

The user assumes full responsibility in the case of operations carried out without observing the instructions in the manual. Primary general safety regulations:

- do not touch the unit with wet hand and/or feet
- do not insert screwdrivers or other pointed objects between guards or moving parts of the appliance
- do not pull the power cord to disconnect the appliance from the electrical mains
- make sure that the appliance is not used by unsuitably qualified persons
- before performing any clearing or maintenance on the appliance disconnect it from the electrical mains by switching of the main switch and extracting the plug
- in the case if faults or malfunctions, switch off the appliance and do not attempt to repair it yourself. All the service and repair operations must be performed exclusively by suitably qualified authorized technicians.

### 3.2 APPLICATIONS AND INTENDED USE

This appliance has been conceived to be used in hospitals, laboratories, pharmacies, etc. on professional premises. The appliance has been designed for the storage of products at a controlled temperature, within the following temperature ranges, according to appliance serie:

MEDIKA	+2°C +15°C
LABOR	0°C +15°C
COMBI	C+ : 0°C +10°C / -3°C +6°C , C- : 0°C +10°C / -5°C -20°C
EMOTECA	+4°C
MINIFRIGO	+2°C +15°C
CHROMATOGRAPHY	+2°C +15°C
VISION	-10°C -20°C
FREEZER	-10°C -20°C
SUPERFREEZER	-20°C -30°C
ARTIC	-15°C -30°C
SUPER- ARTIC	-25°C -40°C
TER	+2°C +30°C
SPARK-PROOF	+2°C +15°C

All the above listed series of appliance are suitable for the products storage, because of that is suggestible to store only products that are already cooled, or frozen, depending from the serie. On request, may be supplied models having temperature ranges different from the standards specified hereabove, or with improved performances, like for the models with tropicalized cooling unit.

All uses except authorized uses of the appliance shall be considered as “improper use” for which the manufacturer declines all responsibility.

### 3.3 USE CONTRAINDICATION

The appliance must not be used:

- Exposing it to outdoors
- With reductions or multi-way adapters
- In places subject to explosive atmosphere or with risk of fire
- Near to heat sources

### 3.4 SAFETY AND ACCIDENT PREVENTION

The appliance embodies various features designed to assure the safety and protect the health of the user. The following list describes the protections adopted mechanical risks:

- **Stability:** the appliance is designed and built so that even with the shelves/drawers fully extracted in the intended conditions of operation, it will remain stable so that it can be used with no risk of tipping, falling or sudden movement.
- **Surfaces, edges, corners:** accessible parts of the appliance have no sharp corners, sharp edges or rough surfaces that could cause injury.
- **Moving parts:** moving parts of the unit are designed, built and configured to avoid risk. Moving parts are protected by fixed guards to prevent accidental contact that could result in injury.

HEREBELOW ARE LISTED MEASURES ADOPTED FOR THE PROTECTION AGAINST ADDITIONAL RISKS:

- o **Electrical power:** the appliance is designed, built and fitted up with the aim of preventing the risk of electric shock in compliance with established safety regulations
- o **noise:** the appliance is designed and built to reduce risk related to the emission of airborne noise to a minimum (inferior than 70 db).

### 3.5 SAFETY PROTECTIONS ADOPTED

**It is strictly forbidden:**

- to tamper with or remove the evaporator cover that protects the user from the risk of cutting on the heat exchanger fins
- to remove the dataplate fixed to the inside edge of the motor housing showing technical specifications and earth connection warning
- to remove the dataplates on the evaporator unit cover near the electrical wiring inside the motor housing which warn the user to disconnect electrical power before working on appliance

**The manufacturer declines all responsibility for safety of the appliance if the above recommendations are not observed.**

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The appliance must be transported and handled exclusively in a vertical position, in observance of the instructions printed on the packing.

This precaution is necessary to avoid contamination of the refrigerant circuit with compressor lube oil with resulting valve and heat exchanger coil failure and problems starting the electric motor. The manufacturer accepts no responsibility for problems due to transport executed in conditions other than those specified above.

The internal equipment (shelf runners, supports, wire shelves, drawers) are shipped inside the unit. The appliance is secured to a wooden base by means of screws and wrapped in polyethylene or packed in a carton, case or crate.

The appliance must be handled using a fork lift truck with suitable forks (forks length at least equal to 2/3 length of unit).

In case the appliance must be laid down flat in order to bring it into the usage place, it is suggestible to wait at least 6 hours before switch the appliance on.

## 4.2 POSITIONING

Incorrect positioning can cause damage to the appliance and generate hazardous conditions for personnel. The installer must therefore observe the following general regulations:

- make sure you maintain a minimum of 3 cm from the walls
- the room must be well ventilated
- keep well away from sources of heat
- avoid direct sunlight
- remove packing material
- remove accessories from the inside the unit
- remove the wooden base
- position the appliance with the help of a spirit level. Adjust the levelling feet on the metal base of the unit if necessary (in the models fitted with adjustable feet) (Fig.1)
- remove the protective PVC film from the external surfaces of the unit

The appliance is thoroughly cleaned in our factory before delivery, in any case, we recommend, before use, to clean the interior with a soft damp cloth and alcohol in order to eliminate the protective oils.


Fig.1

## 4.3 WIRING AND ELECTRICAL HOOK-UP

The electrical plant and electrical hook-up operations must be performed by a qualified electrician. For safety reasons adhere to the following indications:

- check that the electrical plant is suitably sized for the absorbed power of the unit
- if the electrical socket and the plug on the appliance power cord are incompatible, change the plug with a suitable component, ensuring the replacement part is of the approved type
- do not use reductions or multi-way adapters

It is important to connect the appliance correctly to an efficient earth system executed in compliance with the relevant legislation.

## 4.4 SET-UP OPERATIONS

To avoid errors and accidents, perform a series of checks for possible damage sustained during transport, installation and hook-up operations before starting up the unit.

### 4.4.1 Preliminary checks:

- check the condition of the power cord (no cuts or chaffing)
- check that the feet, door hinges and shelf supports are stable
- check the condition of internal and external components (pipelines, heat exchanger elements, fans, electrical components, etc.); check also that all parts are firmly fixed into position
- check that the door seals are not damaged (broken or scratched) and that the doors close and seal properly.

The user must also observe the following instructions to obtain the best operation from the appliance.

#### 4.4.2 Indication for optimal duty:

- do not block the motor compartment air vents
- arrange the material on suitable shelves or drawers. Do not place the products to refrigerate directly on the base or against the walls, doors or fixed guards of the unit
- make sure door is kept closed
- keep the defrost water drain outlet clear
- limit the frequency and duration of opening; each time the door is opened the internal temperature will alter and with the possibility of ice formations on the evaporator
- perform routine maintenance regularly (see condenser cleaning Par. 6.3 )
- load the material at ambient temperature gradually to allow correct refrigeration
- the power supply must be 230V  $\pm$  10%
- The appliance is designed and built to work in ambient temperatures between +15°C e +25°C and relative humidity of 65%. In ambient working operations different respect the specified, will not be possible to achieve the performances declared from the manufacturer.

#### 4.4.3 Special warnings for models Artic e Super Artic

These appliances have been designed and built to store and maintain products (so, it is suggestible to store products already frozen) with a controlled temperature, that at its maximum is -40°C for the serie SUPERARTIC and -30°C for the serie ARTIC and the rotation of the stored products can not exceed the 5% daily.


#### 4.4.4 INTERNAL FITTING NEW MODELS:

<b>MEDIKA</b>	<b>140-170 - 200 -250- 400 - 500-700Lux</b>
<b>MEDIKA</b>	<b>280 2T - 400 2T - 500 2T - 600 2T</b>
<b>LABOR</b>	<b>140-170 - 200 -250- 400 - 500-700Lux</b>
<b>EMOTECA</b>	<b>140 - 170 - 250 - 400 -700Lux</b>
<b>FREEZER</b>	<b>400</b>
<b>ARTIC</b>	<b>400</b>
<b>VISION</b>	<b>400</b>
<b>TER</b>	<b>200</b>
<b>ANTI-SCINTILLE</b>	<b>400</b>


The innovative system fitted of stainless steel racks (standard fitted) allow the possibility of having an internal fitting mixing shelves and drawers (mounted on telescopic slides) perfectly interchangeable. (For eventual orders of further storage shelves/drawers, please refer always to the model and serial number, see data plate paragraph 10.1).

Position the shelves support on the rack at the desired position, insert them into the special hole (Fig. 4) and turn them of 90° to block them. Having positioned 4 supports on the same height level, is possible to insert the shelf (Fig.3).

In order to change the position of the drawers is needed to extract the drawer itself, and once it is fully extracted, unlock it from the slides through the unlocking devices (plastic lever black colour) positioned at the side. The unlocking of the drawer can be done pushing the right lever up and the left lever down at the same time. After that is possible to remove the telescopic slides lifting up the front part (Fig.2) of the slide itself, in order to extract the slide from the rack in the front part and pulling the slide forward, in order to extract the rear part. Repeat the same operations, but inverting the sequence to reposition the slides and the drawers.


**Fig.2**

**Fig.4**

**Fig.3**

The wall-racks are even easily removable in order to allow a better cleaning of the interior of the appliance. The racks mounted in the front part can be removed pushing them upward (Fig.5), while the ones mounted in the rear side can be removed first, unscrewing safety fixing screw (using a cross screwdriver) positioned in the top part (Fig.6), and then pushing against the internal wall the small plastic edge on the top part of the rack and simultaneously pulling the rack itself upward (Fig 5).

**Fig.5**

**Fig.6**

#### 4.4.5 INTERNAL FITTING MODELS:

MEDIKA	600 – 700 – 900 – 1000 – 1500 – 2300
LABOR	600 – 700 – 900 – 1000 – 1500 – 2300
FREEZER	600 – 700 – 900 – 1000 – 1500 – 2300
SUPERFREEZER	700
ARTIC	600 – 700 – 1000 – 1500 – 2300
SUPER-ARTIC	700
SUPER-ARTIC	700 2T
COMBI	700 2T – 1000 2T – 1500 2T – 2300 2T
CROMATOGRAFIA	700
ANTI-SCINTILLE	600 – 700 – 1000 – 1500
EMOTECA	700 – 1500
TER	700

For these models the mounted on the shelves is different respect what described in the previous paragraph.

To fix the shelves correctly, position the shelves slides at the desired position, insert first the slide into the wall-rack mounted on the internal rear side and then insert the edge of the side of the slide into the wall-rack mounted on the internal front part at the corresponding height of the fixing in the rear. Repeat the same operation for both the slides paying attention that are not interchangeable, there are right slides and a left slides. Finally, having inserted a pair of slides at the same height, insert the shelf between them (Fig.7).



Fig.7

#### 4.4.6 Internal fitting with brackets

For the models Medika 130, 270, 450 2t , Ter 80, Labor/Medika/Freezer 900 the shelves are supported by brackets.

### 5.1 COMPLETE CONTROL PANEL

	1	TEMPERATURE CHART RECORDER : recording continuously the internal temperature of the appliance
	2	ALARM SYSTEM WARNING LAMPS (RED WARNING LAMPS) : advising in case the internal temperature is above or below the programmed working temperature
	3	ALARM SYSTEM WARNING LAMP (ORANGE WARNING LAMPS) : advising in case of power failure
	4	ALARM SYSTEM SWITCH (ORANGE SWITCH) : test button of the alarm system
	5	ALARM SYSTEM MAIN SWITCH (GREEN SWITCH) : to switch the alarm system ON
	6	THERMOSTAT : allowing the programming of the working temperature and displaying the actual internal temperature of the appliance
	7	LIGHTING SWITCH (GREEN SWITCH) : to switch on the internal lighting, <u>OR</u> HUMIDITY SWITCH (BLUE SWITCH) : to adjust the internal humidity level, H=high, L=low, <u>OR</u> DEFROSTING LAMP (ORANGE PILOT LAMP) : advising when the defrosting is in progress, <u>OR</u> BLACK COVER : pre-arrangement.
	8	MAIN SWITCH (GREEN SWITCH) : to switch the appliance ON

**NOTE:** The **ALARM SYSTEM** (points nr. 2,3,4,5) and the **TEMPERATURE CHART RECORDER** (point nr. 1) are optionally available, so, if not ordered they may not be present

The appliances may be provided, on request, with the main switch (point nr. 8) key operating.

## 5.2 THERMOSTAT

The refrigerator is fitted up with an electronic and digital (mod. ID 971 or ID974 or IS 972) depending from the model.

In case the appliance is fitted with a different thermostat respect the ones mentioned above, the specific instruction manual is provided

**ATTENTION** : For the models with 2 temperatures with the doors one above the other: Models **MEDIKA 400 2T** and **MEDIKA 500 2T** the right hand side thermostat controls the bottom compartment, while the left hand side thermostat controls the top compartment.

Models **LABOR 600 2T – LABOR 700 2T** the right hand side thermostat controls the top compartment, while the left hand side thermostat controls the bottom compartment.

The user has a display and four keys for controlling status and programming of the instrument:

<b>UP KEY</b>		Scrolls through the menu items Increases the values Activates manual def. function
<b>DOWN KEY</b>		Scrolls through the menu items Decreases the values Programmable by parameter
<b>FNC KEY</b>		ESC function (exit) Programmable by parameter
<b>SET KEY</b>		Accesses the set point Accesses the menus Confirms the commands Displaying the alarms (if actives)


At start-up the instrument performs a Lamp Test; for few seconds the display and the leds blink, in order to verify their integrity and correct operation.

### Temperature setting

Access the “Machine Status” menu by pressing and quickly releasing the “set” key. The label of the “set” folder appears. To display the Set point value press the “set” key again. The value appears on the display.

To change the Set point value, use the “UP” and “DOWN” keys within 15 seconds.

### Temperature alarm signalling (alarm incorporated among the thermostat functions)

The alarm condition is always signalled by the buzzer and by the led of alarm icon ().

The regulation of the maximum and minimum temperature alarm refers to the room probe. The temperature limits are defined by the parameters “HAL” (maximum alarm) and “LAL” (minimum alarm). The alarm limits are factory setted, in case they have to be modified or further information/clarification would be necessary, please refer to the attached thermostat booklet.

THIS FUNCTION IS DISABLED IN CASE THE APPLIANCE IS EVEN PROVIDED WITH THE VISUAL AND ACOUSTIC ALARM SYSTEM.

### Other alarm signalling (alarm incorporated among the thermostat functions)


The alarm signal produced by a faulty room probe (probe 1) is shown as “E1” on the instrument display.

The alarm signal produced by a faulty evaporator (probe 2) is shown as “E2” on the instrument display.

### Manual activation of the defrosting cycle


If might be required, it is possible to manually activate the defrosting cycle, pressing the “UP” key for 5 seconds. If defrosting conditions are not present, (for example the evaporator probe temperature is higher than defrost stop temperature), the display will blink three (3) times, in order to indicate that the operation will not be performed.

**A small LED lights up, under the symbol when the function is activated.**

	Compressor or relay 1	ON for compressor ON, flashing for delay, disabled protection or activation
	Defrost	ON for current defrosting; flashing for manual activation
	Alarm	ON for an active alarm; flashing for silenced alarm
	Fans	ON for an active fan

**CAUTION: NEVER ATTEMPT TO MODIFY THE PROGRAMMED WORKING PARAMETER OF THE THERMOSTAT AS THE CORRECT REFRIGERATOR WORKING CAN BE COMPROMISED.**

### **5.3 VISUAL AND ACOUSTIC ALARM SYSTEM (IF PRESENT)**

The appliance can be provided with a visual and acoustic alarm system for temperature deviations and power failure, the system is battery operating in order to grant the alarm signalling even in case of lack of mains and the battery is kept continuously charged during the normal working conditions (presence of the mains). The alarm system is compounded by switches and warning lamps fitted on the control panel (points nr. 2,3,4,5 of the paragraph “Control Panel”) and from the alarm system box (see above) that is located in the motor compartment (in the top or in the bottom part of the cabinet, depending from the model). The alarm system box is also provided with 2 outputs for the remoting of the alarm signal, one is countersigned with a label “12V” as provide a 12 Volts contact, while the second is countersigned with a label “NO/NC” as provide a potential-free contact.

Before to switch ON the alarm system main switch (point nr. 5 of the paragraph “Control Panel”), please wait until that the refrigerator has achieved the temperature programmed

through the thermostat (see paragraph 'Thermostat', point nr. 6 of the paragraph "Control Panel").

At start-up the instrument performs a lamp and a buzzer test; for few seconds the display of the alarm system box blinks (showing EEE) and the buzzer sounds in order to verify the integrity and correct operation of the alarm system itself.

The alarm system setpoints (both the minimum and maximum values) are factory setted considering a working temperature for the appliance, but, in any case they remain completely adjustable according to the user needs. The adjustment of the alarm system setpoints can be done using the alarm system box: pressing the 'SET' button and it will appear on the display 'HI' (high temperature alarm, i.e. +10°C) or 'LO' (low temperature alarm, i.e. 0°C) then program the temperature using the buttons '+' (to increase the value) or '-' (to decrease the value) until reaching the desired value. The new value will be automatically memorized after 5 seconds from the last change.

After 60 seconds from the last change, the display will enter in "energy saving" condition (to save the battery electricity) and it displays '.' (dot).

Please note that is absolutely essential to keep a minimum of 3°C difference for both the high temperature alarm and low temperature alarm the respect the setpoint temperature programmed on the thermostat (see paragraph 'Thermostat', point nr. 6 of the paragraph "Control Panel"), in order to avoid alarm signalling due to minimum variation of the internal air temperature that do not effect the stored products (i.e. door openings). Example for alarm settings:

FOR A POSITIVE TEMPERATURE CABINET (ABOVE 0°C):

SET POINT	HI	LO
0°C	+ 5°C	-2°C
+ 1°C	+ 6°C	-1°C
+ 2°C	+ 7°C	0°C
+ 3°C	+ 8°C	0°C
+ 4°C	+ 9°C	+1°C
+ 5°C	+10°C	+2°C
+ 6°C	+10°C	+3°C
+ 7°C	+11°C	+4°C
+ 8°C	+12°C	+6°C
+ 9°C	+13°C	+7°C
+10°C	+11°C	+8°C
+11°C	+15°C	+9°C
+12°C	+16°C	+10°C
+13°C	+17°C	+11°C
+14°C	+18°C	+12°C
+15°C	+19 C	+13°C

FOR A NEGATIVE TEMPERATURE CABINET (BELOW 0°C):

SET POINT	HI	LO
- 1°C	+ 6°C	- 3°C
- 10°C	+ 0°C	- 12°C
- 15°C	- 5°C	- 18°C
- 20°C	- 10°C	- 23°C
- 25°C	- 15°C	- 28°C
- 30°C	- 20°C	- 34°C
- 40°C	- 30°C	- 43°C

It is possible to test even the alarm system integrity and correct operation (i.e. that the battery is charged, the buzzer works, etc.) pressing the test switch on the control panel (point nr. 4 of the paragraph “Control Panel”).

The following alarm conditions are signalled:

- Appliance internal temperature higher than the programmed : light-up the symbol  (point nr. 2 of the “Control Panel”) and the buzzer sounds
- Refrigerator internal temperature lower than the programmed : light-up the symbol  (point nr. 2 of the “Control Panel”) and the buzzer sounds
- Power failure : light-up the orange colour warning lamp  (point nr. 3 of the “Control Panel”) and the buzzer sounds

#### Alarm system sensor failure signalling

In case of failure of the alarm system sensor, the display will show ‘EEE’, the acoustic alarm will activate and the warning lamps on the front panel (point nr. 2 of the paragraph “Control Panel”) will blink alternatively. In this case proceeds replacing the alarm system/thermostat sensor (both the thermostat and alarm system use the same sensor).

## **5.4 TEMPERATURE CHART RECORDER (IF PRESENT)**

The appliance can be provided with a temperature chart recorder with 7 days recording continuously the internal temperature of the appliance itself on special round charts (125 mm Ø). The charts have a temperature range allowing the recording of all the possible temperature deviations and they are supplied in a quantity of 100 pieces to permit a recording period near to the 2 years.

The temperature on the charts is lead by a special ink cartridge, with blue colour ink.

To prevent undesired tampering, the chart is protected by a transparent plexiglass door, provided with key locking.

The temperature chart recorder is fed by a 1,5V battery in order to grant the recording even in case of power failure.

### **5.4.1 Changing of the paper chart (Fig.11)**

Release the pointer from the diagram, unscrew the knurled nut (4), place the diagram on the drive’s axis, paying attention to insert the diagram (3) under the tabs foreseen to hold it, and tighten the nut, the knurled part toward the outside, up to the stop, finally, gently bring the pointer on the diagram on the exact day of recording.

### **5.4.2 Replacing of the battery(Fig.12) (Battery LR6 / AA 1,5 V)**

Open the recorder door, lift the pointer so as to release the diagram, unscrew the knurled nut and remove the diagram. Firmly hold the case in one hand and, with the other, pull the mechanism axis toward the outside by alternatively swinging up and down to release it from its base. Change the battery located behind the mechanism. The mechanism thus taken out allows access to the speed-changing lever, with the possibility to select the daily recording (1 x 24h) or the weekly recording (7 x 24 h). Repeat the same operations, but inverting the sequence to reposition all the components.


Fig.11



Fig.12

### 5.4.3 Replacing of the ink-tip

Lift the plate's pointer. Insert the end of the ink-tip into the fibre point pen's slide rail up to the stop. Gently place the pointer on the diagram again

**IMPORTANT: Do not place the ink-tip in contact with the fingers**

## 6

## ORDINARY CARE AND MAINTENANCE

The informations in this section regard the end user, or other non-specialized personnel, and routine maintenance technician.

### 6.1 PROHIBITION OF THE REMOVAL OF THE SAFETY DEVICES

It is strictly forbidden to remove guards or safety devices when performing routine maintenance works. The manufacturer disclaims all liability that may arise if this regulation is not observed.

### 6.2 CLEANING THE INTERIOR AND EXTERIOR OF THE APPLIANCE

The appliance is thoroughly cleaned in our factory before delivery. We recommend, however, that you clean the interior of the appliance before use. Before any cleaning operation, make sure that the appliance power cord is disconnected. Also it is suggestible to clean both the interior and the exterior surfaces of the appliance at least twice in a year.

For more details see the below paragraph.

For this scope is suggestible :

- Cleaning products: water and non-abrasive neutral detergent. **DO NOT USE SOLVENT OR THINNERS**
- Cleaning method: use a cloth or sponge soaked in a suitable cleaning product to clean the interior and exterior parts of the cabinet

- Sanitation: do not use substances that could alter the basic characteristics of the stored material
- Rinsing: use a cloth or sponge soaked in clean water. **DO NOT USE WATER JETS**
- Frequency: at least twice in a year or at different intervals in accordance with the type of pharmaceutical products conserved

### 6.3 CONDENSER CLEANING


In the models with the motor in the bottom part remove the protection guards as described in the paragraph 7.1.1

In the models with the motor in the top part, the condenser is directly accessible using a step ladder to reach the top of the appliance.

Clean with MONTHLY frequency the condenser (depends from the dust present in the ambient). Remove any dust or fluff that has deposited on the heat exchanger fins using an air jet or dry brush and, working with up and down movements. **DO NOT USE METAL BRUSH**

#### **ATTENTION:**

**BEFORE CLEANING THE CONDENSER SWITCH OFF THE APPLIANCE, DISCONNECT THE POWER CORD.**


**In order to guarantee an optimal working of the appliance is indispensable follow the indication given from the manufacturer and arrange the ordinary maintenance through qualified technician.**

### 6.4 CONDENSATE WATER DRAINING

The appliance is designed and built with the arrangement of automatic (cyclic) defrostings, factory pre-programmed through the thermostat. The defrostings originate water.

In the models with the motor in the top part, the water is drained into a tray, including into the standard supplying, that has to be positioned under the appliance itself, supported by pre-arranged slides.

This tray has to be emptied cyclically, the frequency depends from the water drained.


List of the models with the motor in the top part:

MEDIKA	600 – 700 – 900 – 1000 – 1500 – 2300
MEDIKA	450 2T – 1000 2T – 1500 2T – 2300 2T
LABOR	400 – 450 – 600 – 700 – 900 – 1000 – 1500 – 2300
FREEZER	400 – 450 – 600 – 700 – 900 – 1000 – 1500 – 2300
ARTIC	400 – 450 – 600 – 700 – 1000 – 1500 – 2300
SUPER-ARTIC	700
SUPER-ARTIC	700 2T
COMBI	600 2T – 700 2T – 1000 2T – 1500 2T – 2300 2T
CROMATOGRAFIA	700
ANTI-SCINTILLE	400 – 600 – 700 – 1000 – 1500 – 2300
EMOTECA	700 – 1500 – 2300

Optionally available, for the models with the motor in the top part, there is the possibility to have a special heated tray (electrically operating) that evaporates the defrosting water automatically.

In all the models with the motor in the bottom part evaporate the defrosting water automatically.

Also, there are some special versions, having the cooling unit “monoblock type”, located in the top part, that evaporate the defrosting water automatically, without the need of any extra device.

## 7

## EXTRAORDINARY MAINTENANCE AND REPARATIONS

Extraordinary maintenance and reparations are to be performed exclusively by specialized technicians authorized by the manufacturer.


The manufacturer declines all liability in the case of works performed by the user or unauthorized persons, or if non-original spare parts are fitted to the appliance.

All malfunctioning due to the fitting of non-original spare parts will not be accepted from our technicians and will cause the immediate loss of the **GUARANTEE**.


### 7.1 PROTECTIONS REMOVAL (for the models the motor in the bottom part).

To perform extraordinary maintenance and reparations is needed to remove correctly the protection guards as described below:


### 7.1.1 Bottom guard

**Fig.13**

**Fig.14**

**Fig.15**

**Phase 1** Open the appliance door e remove, using a cross screwdriver, the 3 screws fixing the bottom guard.(Fig.13)

**Phase 2** Close the door. The guard is blocked in the bottom part with 2 “hooks”, so, must be inclined frontward (Fig.14) and slide it out from the top(Fig.15).

## 7.1.2 Top front panel and electric circuit cover


Fig.16


Fig.17


Fig.18

To remove the control panel and have the access to the electric circuit follow the below listed operations:

**Phase 1** Open the appliance door and remove, using a cross screwdriver, the 3 screws fixing the control panel.

**Phase 2** The control panel blocked in the bottom part with 2 “hooks”, so, must be inclined frontward (Fig.16) and slide it out to the right side (Fig.17).

**Phase 3** To remove the top cover, unscrew, using a cross screwdriver, the fixing screws of the rear-side guard and push the top cover frontward until its complete unblocking from the fixing hooks (Fig.18).

Repeat the same operations, but inverting the sequence to reposition all the components.

8

DEMOLITION

The product has been realised with recyclable material.  
This appliance complies with the 2002/96/EC European Directive.

The symbol  on the product means that it must not be considered as a domestic waste but it must be handed near the competent authority that recycles electric and electronic appliances.

Before scrapping the machine, make it unusable. First of all cut the connecting cable, remove the doors, tiers and drawers to do not let children go inside the product.

Do not waive it even for a few days.

For further information about management, retrieval and recycle of the product, please contact the local office, the domestic wastes picking service or the distributor.

Respect the laws applicable.

The refrigerant gas present into the cooling circuit must be extracted by authorized personnel.

**9**

**ATTACHMENTS**


The following documents are attached:

- Declaration of conformity with the directive 73/23 CEE
- Declaration of conformity with the directive 89/336 CEE
- Declaration of conformity with the directive 89/109 CEE
- Declaration of conformity to the standard CEI EN 61010-1 (CEI 66/5)
- Production and test sheets
- Electrical circuit diagrams


**10**

**DATA PLATES**

**10.1 DATA PLATES WITH THE CHARACTERISTICS OF THE APPLIANCE**

**10.2 OTHER INDICATION LABELS**

 <p><b>PRIMA DI APRIRE LA PROTEZIONE TOGLIERE LA TENSIONE</b>  <b>TURN OFF AND UNPLUG AC BEFORE OPENING COVER</b>  <b>AVANT D'OUVRIER LA PROTECTION ÔTER LA TENSION</b>  <b>BEVOR DER SCHUTZ ZU OEFFNEN, ZU ENTSPANNEN</b></p>	<p>Turn off and and unplug AC before opening covers</p>
	<p>Periodic condenser cleaning</p>
	<p>Earthing</p>

In the below sheet is listed the most common faults, the possible causes and the relative action to be taken.

FAULT	POSSIBLE CAUSE	ACTION
Appliance does not switch on	Switch set to 'OFF' power failure  Other	Set switch to 'ON' Check plug, socket, fuses, electrical line Contact service department
Refrigeration unit does not start	Defrosting cycle in progress  Control panel fault Other	Wait for cycle end, switch off and then on again Contact service department Contact service department
Refrigeration unit does not cut out at set temperature	Control panel fault Temperature sensor defective	Contact service department Contact service department
Refrigeration unit runs constantly but does not reach set temperature	Room temperature too high Condenser dirty Evaporator is coated with ice Refrigerant needs recharging Condenser fan not running Evaporator fan not running Door seals are not efficient  Defrost valve is open	Provide better ventilation of room Clean condenser Defrost manually Contact service department Contact service department Contact service department Check seals / make sure refrigerator contents not obstructing doors Contact service department
The refrigeration unit is noisy	Appliance not levelled  Moving parts touching guards  Screws unscrewed Other	Check that the appliance is levelled Check that the fans do not touch their guards or pipes Fix the screws Contact service department
Water or ice deposits in the drip tray	Drain outlet is clogged Appliance not level	Clean drain and drain outlet Contact service department
Build up of ice on evaporator	Improper use Control panel fault Door seals are not efficient	Contact service department Contact service department Check seals / make sure refrigerator contents not obstructing doors
The display of the alarm system box is switched off	Alarm system main switch set to 'OFF' Fuse blown Battery discharged Battery failure Other	Set switch to 'ON' Sontact service department Contact service department Contact service department Contact service department

**In order to be able to arrange a prompt intervention, at the moment of the calling, please provide the model of the appliance the the relative serial number that can be found on the data plate sticked in the rear side of the appliance or on the present manual**

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