Bartels Micropumps


Micropumps transporting the tiniest amounts of gases or liquids can be considered the heart of microfluidics.

In many sectors they have become indispensable. Dosing lubricants, feeding sensors with sample gas or mixing starch into the steam of flat irons are only a few of the manifold tasks they can fulfill. Many further fields of application for example are located in medical technologies and analytics.


Extremely small in size and low in weight, with good particle tolerance and temperature resistance, Bartels micropumps are well prepared to be used in any of these sectors. As they are almost completely made of plastics, large quantities of these pumps can be produced at low cost and so may well be used as disposables.

The functional principle of the Bartels micropumps is based on a piezoelectric diaphragm in combination with passive check valves. A piezo ceramic mounted on a coated brass membrane is deformed when voltage is applied. By the resulting down stroke, the medium is being displaced out of the pump chamber below. The check valves on both sides of the pump chamber define the flow direction. When the voltage decreases, the corresponding piezo deformation causes an upstroke of the membrane. The medium is sucked in and the chamber is filled again. In every second, the pump can do several hundreds of such pumping cycles. The pumping performance can be influenced by adjustment of the parameters.


Important advantages for all users result from the radically simple pump design: Injection molded parts for housing and pump chamber, piezo actuators and passive valves constitute the key components. Thus any adaptation to specific requirements concerning flow rate or back pressure is easy to realize. This customization of micropumps with the appropriate electronic controllers is part of the services offered by Bartels microComponents. If requested, the pumps can be fully integrated into complex system designs as well.

mp6 – double actuator micropump

mp5 – smallest available plastic micropump

Functional principle.


Once the perfect pump for your application has been found, you may purchase an exclusive production license for this version to include the component into your own production processes. Of course Bartels microComponents can also realize a high quality serial production for you at low cost.

mp6 Micropump

The Bartels micropump mp6 combines two piezo actuators inside a single housing. The new pump joins the established functional principle and central advantages of its parent generation mp5 with its own **specific innovative features**.

The small power pack can handle **twice the back pressure** the mp5 can cope with, has an **increased priming capability** and is of **higher bubble tolerance**, so that even gas-liquid-mixtures can be pumped without problems. Its **low power consumption** is a further advantage. In the entire pump only one material gets into contact with the medium. Now two types of materials are available. In the standard mp6 all relevant parts are made of polyphenylsulfone (PPSU) and in the mp6-pp of polypropylene (PP).


The mp6 can already be offered at low prices for large quantities due to an automated assembly. The mp6-pp is produced less automated but for larger quantities it can also be transferred into serial production.

mp5 Micropump


The mp5 from Bartels microComponents is the smallest and lightest micropump available. Since 2004 the mp5 has successfully shown the potential of piezo membrane pumps. Due to the limited bubble tolerance and higher price in medium quantities, it is now replaced by its successor the mp6 in many applications. If size is a challenging factor, then the mp5 is still the right choice.

Its low power consumption and tiny size makes the mp5 the perfect pump to be fully integrated into your product's design. Test the mp5 now and ask about the possibilities of a customer specific adaptation - for your individual micropump.


Typical characteristics of the Bartels micropumps:


volume flow versus pressure

volume flow versus frequency

volume flow versus amplitude

Technical Data of the mp6¹

mp6	Order code: mp6
Pump type	piezoelectric diaphragm pump
Number of actuators	2
Dimensions without connectors	30 x 15 x 3.8 mm 1.1811 x 0.5906 x 0.1498 in.
Weight	2 g
Fluidic connectors	barbed tube clip, (outer diameter 1.9 mm, length 3.5 mm) ²
Electric connector	flex connector 1.25 mm pitch
Power consumption	< 200 mW
Self-priming	yes ³
Pumping media	liquids, gases, mixtures
Operating temperature	0–70°C
Life time	5000 h ⁴
IP code	IP33 ⁵
Material in contact with media	polyphenylsulfone (PPSU) ⁶
Suitable evaluation controller	mp-x, mp6-EVA and mp6-OEM
Typical values of flow and back pressure for selected media (values defined with mp-x: 250 V, SRS):	
Gases (as product mp6-air)	
max. volume flow	18 ml/min (300 Hz) ⁷
max. back pressure	100 mbar (300 Hz) ⁷
Liquids – water	
max. volume flow	7 ml/min +/- 15% (5.95 ml/min – 8.05 ml/min) ⁴
max. back pressure	600 mbar +/- 15% (510 – 690 mbar) ⁴

¹ Typical values. Values can vary under application conditions. Content is subject to changes without notice.

² Recommended tubing: 1.3 mm inner diameter.


³ Conditions: Suction pressure < 10 mbar, DI water, settings mp-x: 100 Hz, 250 V, SRS, the max. volume flow will be reached after a few minutes of operation time.

⁴ Conditions: DI water, room temperature, settings mp-x: 100 Hz, 250 V, SRS

⁵ Can be changed to IP44.

⁶ The mp6 is not stable against concentrated alcoholic solutions as MeOH or EtOH.

⁷ Conditions: Gases, room temperature, mp-x: 300 Hz, 250 V, SRS



Please find more information concerning the controller and the equipment in the corresponding data sheets.

Technical Data of the mp6-pp¹

mp6-pp	Order code: mp6-pp
Pump type	piezoelectric diaphragm pump
Number of actuators	2
Dimensions without connectors	30 x 15 x 3.8 mm 1.1811 x 0.5906 x 0.1498 in.
Weight	2 g
Fluidic connectors	barbed tube clip, (outer diameter 1.9 mm, length 3.5 mm) ²
Electric connector	flex connector 1.25 mm pitch
Power consumption	< 200 mW
Self-priming	yes ³
Pumping media	liquids, gases and mixtures
Operating temperature	t.b.d.
Life time	t.b.d.
IP code	IP 33 ⁴
Material in contact with media	polypropylene (PP) ⁶
Suitable evaluation controller	mp-x, mp6-EVA and mp6-OEM
Typical values of flow and back pressure for selected media (values defined with mp-x: 250 V, SRS):	
Gases	
max. volume flow	14 ml/min +/- 15% (11.9 – 16.1 ml/min) ⁵
max. back pressure	150 mbar +/- 15% (127.5 – 172.5 mbar) ⁵
Liquids – water	
max. volume flow	5 ml/min +/- 15% (4.25 – 5.75 ml/min) ⁶
max. back pressure	650 mbar +/- 15% (552.5 mbar – 747.5 mbar) ⁶
Liquids – MeOH	
max. volume flow	6.8 ml/min ⁷
max. back pressure	550 mbar ⁷

¹ Typical values. Values can vary under application conditions. Content is subject to changes without notice.

² Recommended tubing: 1.02 mm inner diameter.

³ Conditions: Suction pressure < 10 mbar, DI water, settings mp-x: 100 Hz, 250 V, SRS, the max. volume flow will be reached after a few minutes of operation time.

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⁵ Conditions: Gases, room temperature, mp-x: 300 Hz, 250 V, SRS

⁶ Conditions: DI water, room temperature, settings mp-x: 100 Hz, 250 V, SRS



⁷ Conditions: DI water, room temperature, settings mp-x: 100 Hz, 250 V, SRS

Please find more information concerning the controller and the equipment in the corresponding data sheets.

Technical Data of the mp5¹

mp5	Order code: mp5
Pump type	piezoelectric diaphragm pump
Number of actuators	1
Dimensions without connectors	14 x 14 x 3.5 mm 0.5512 x 0.5512 x 0.1378 in.
Weight	0.8 g
Fluidic connectors	barbed tube clip, (outer diameter 2 mm, length 3 mm) ²
Electric connector	flex connector / phone jack
Power consumption	< 200 mW
Self-priming	yes ³
Pumping media	liquids or gases
Operating temperature	0 – 70°C
Life time	5000 h ⁴
IP code	IP44
Material in contact with media	polyphenylsulfone (PPSU), polyimide (PI), nitrile butadiene rubber (NBR) ⁵
Suitable evaluation controller	mp-x, mp6-EVA and mp6-OEM
Typical values of flow and back pressure for selected media (values defined with mp-x: 250 V, SRS):	
Gases	
max. volume flow	15 ml/min (300 Hz) ⁶
max. back pressure	30 mbar (300 Hz) ⁶
Liquids – water	
max. volume flow	5 ml/min (100 Hz) ⁴
max. back pressure	250 mbar (100 Hz) ⁴

¹ Typical values. Values can vary under application conditions. Content is subject to changes without notice.

² Recommended tubing: 1.3 mm inner diameter.

³ Conditions: Suction pressure < 10 mbar, DI water, settings mp-x: 100 Hz, 250 V, SRS, the max. volume flow will be reached after a few minutes of operation time.

⁴ Conditions: DI water, room temperature, settings mp-x: 100 Hz, 250 V, SRS

⁵ The mp6 is not stable against concentrated alcoholic solutions as MeOH or EtOH.

⁶ Conditions: Gases, room temperature, mp-x: 300 Hz, 250 V, SRS

Please find more information concerning the controller and the equipment in the corresponding data sheets.


Accessories for the Bartels Micropumps

mp-x controller	Order code: mp-x
Access to the full range of driving parameters. A system for the professional evaluation of the micropumps. <i>(excl. mp6-connection cable, please order when required)</i>	
Dimensions	7.5 x 16 x 20 cm 2.983 x 6.299 x 7.874 in.
Weight	ca. 800 g
Pumping media	liquids, gases
Max. volume flow (DI-water)	mp6: 7 ml/min (SRS, 250 V, 100 Hz) mp6-pp: 5 ml/min (SRS, 250 V, 100 Hz) mp5: 5 ml/min (SRS, 250 V, 100 Hz)
Adjustable parameters	Amplitude, frequency, signal form
Amplitude range	0 – 250 V
Frequency range	0 - 300 Hz
Signal form	SRS, rectangular, sine
Power supply	mains adaptor
Current consumption	750 mA at 7.5 V
USB-Port incl. driver	one; a CD with driver software is included
connectable micropumps	mp5: 1x-2x or mp6 respectively mp6-pp: 1x

mp-x


mp6-OEM controller	Order code: mp6-OEM
The OEM-controller drives the micropump at adjustable performance in a package similar to an integrated circuit. It enables integration into system electronics or on a PCB.	
Dimensions	10.5 x 20.5 x 6 mm 0.4134 x 0.8070 x 0.2362 in.
Pumping media	liquids, gases
Max. volume flow (DI-water)	4.5 ml/min
Adjustable parameters	Amplitude, frequency
Amplitude range	85 – 235 V
Frequency range	25 - 226 Hz (Frequencies up to 1000 Hz are possible; please note that the amplitude will decrease in that case. Frequencies down to 1 Hz are possible with an external clock generator.)
Signal form	similar to rectangular
Power supply	2.5 – 5.5 VDC (5 V recommended for optimized performance)
Pumping media	ca. 30 mA at 5 V
Pin arrangement	DIL 14; horizontal ~2.54 mm, vertical ~7.62 mm

mp6-OEM

mp6-EVA evaluation board Order code: mp6-EVA	
The evaluation board enables the simple use of the micropump based on the mp6-OEM controller. Next to preset standard parameter (235 Vpp, 100 Hz) the mp6-EVA also allows adjusting the pump parameters flexibly, this partly by external tuning. As the supply voltage of the module can be provided via USB (no data interface), just attach it to a USB power supply and start the evaluation. Alternatively it can also be supplied by a 2.5 – 5 V voltage source.	
Dimensions	6.5 x 3 x 2 cm 2.5590 x 1.1811 x 0.7874 in.
Pumping media	liquids, gases
Max. volume flow (mp6, DI-Water)	4.5 ml/min
Adjustable parameters	Amplitude, frequency
Amplitude range	100 – 235 Vpp
Frequency range	25 – 226 Hz (Frequencies up to 1000 Hz are possible; please note that the amplitude will decrease in that case. Frequencies down to 1 Hz are possible with an external clock generator.)
Signal form	similar to rectangular
Power supply	via USB or 2.5 V – 5.5 VDC
Current consumption	ca. 30 mA at 5 V


mp6-EVA

mp6-con connection cable Order code: mp6-con	
Connector for mp6/mp6-pp to mp-x	
Design and connectors	- Molex FCC 1.25 mm pitch - 85 cm (33.465 in.) cable - Binder 620 connector

mp6-con

mp6-mol connector Order code: mp6-mol	
Connector to micropump mp6/mp6-pp for custom made cabling	
Type	Molex FCC 39532045 1.25 mm pitch
packaging unit	10 pieces

mp6-mol

mp-cv check valve Order code: mp-cv	
The passive check valve eliminates the back flow of the pumping medium, when the micropump is switched off. It can be connected via tubing.	
Dimensions	21 mm x 5.5 mm (length x wrench size) 0.82677 x 0.2165 in.
Materials in contact with the pumped media	Silicone, stainless steel
Fluidic connectors	barbed tube clip, length : 5.6 mm for tubing with internal diameter: 1.3 mm
Cracking pressure	typical < 35 mbar
Max. back pressure	500 mbar
Typical leak rate	<20 µl/h for DI-water (at 500 mbar)


mp-cv

mp-t tubing	
Inlet/outlet compatible Tygon® tubing	
Inner diameter 1.3 mm	Order code: mp-t ID1.3
Inner diameter 1.02 mm	Order code: mp-t ID1.02

mp-t

mp-y tubing connector		Order code: mp-y
Y-connector for tubing, for the parallel use of two micropumps:		
Material	polypropylene (PP)	
for tubing inner diameters of	1.3 – 2.6 mm 0.0512 – 0.1024 in.	

mp-y

Evaluation Sets

mp5-go! Set	Order code: mp5-go!
The evaluation of the mp5 can be started directly with this set.	
It contains:	
3 micropumps mp5	
1 mp-x controller	
1 meter mp-t ID1.3 tubing	

mp6-go! Set	Order code: mp6-go!
The evaluation of the mp6 can be started directly with this set.	
It contains:	
3 micropumps mp6	
1 mp-x controller	
1 mp6-con	
1 meter mp-t ID1.3 tubing	


mp6-go! Set

mp6-basic Set	Order code: mp6-basic
The evaluation of the mp6 can be started directly with this set.	
It contains:	
3 micropumps mp6	
1 mp6-EVA evaluation board	
1 meter mp-t ID1.3 tubing	

mp6-basic Set

mp6-pro Set	Order code: mp6-pro
The evaluation of the mp6 can be started directly on the customer's circuit board with this set.	
It contains:	
5 micropumps mp6	
5 mp6-OEM controller	
1 packaging unit of mp6-mol	
1 meter mp-t ID1.3 tubing	

mp6-pro Set

mp6-pp-go! Set	Order code: mp6-pp-go!
The evaluation of the mp6-pp can be started with this set.	
It contains:	
3 micropumps mp6-pp	
1 mp-x controller	
1 mp6-con	
1 meter mp-t ID1.3 tubing	

The offered accessories and kits are meant to assist your evaluation process. After the feasibility of the micropump in the customer specific application has been proven, an adequate miniaturization of the controller and the equipment can be carried out.

The design of customized controllers is part of the services offered by Bartels microComponents.

Please contact us, if we can support you in choosing the suitable equipment.

Bartels Mikrotechnik Newsletter #9

One for four – Electronic evaluation board for four micropumps

Content:

- mp6-QuadEVA – a new evaluation board to control four micropumps simultaneously and to achieve higher frequencies
- Why more frequency and why more pumps when pumping gas?

mp6-QuadEVA – a new evaluation board to control four micropumps simultaneously and to achieve higher frequencies

As a developer and producer of micropumps we hear about customer applications with our products but also the requests for enhanced performance. One widespread goal of customers over a variety of different application fields is to create small mobile devices. Often a conventional pump – even if it is small for its kind – does not match the available space. Our mp6 micropumps are small, having quite powerful performance but are not always as strong as the customer wants. More performance is quite often a demand for many applications. Sure there are other pumps on the market that can create more flow and/or more pressure, but those are either larger or require much more power and most of them are more expensive.

Though you can easily drive our micropumps with the standard electronics – the lab electronic mp-x and the mp6-OEM for integration in your own circuit – both are not designed to handle more than one pump at the same time. Therefore, if you want to drive two or more pumps in a serial connection to multiply the pressure generation or in parallel connection to raise the flow rate, you will need to have the same number of electronics as you have pumps; one for one.

So, even if the small dimensions of the pumps allow having more than one pump in a mobile device, the need for increasing the amount of electronic drivers may hinder this approach to be reasonable.

Now we have developed a new electronic to address this draw-back: the mp6-QuadEVA.

This evaluation board is not based on the mp6-OEM electronic, as this can only drive one pump at full power. With its own electronic layout the mp6-QuadEVA is meant to control up to four pumps simultaneously. This new approach for controlling multiple pumps was primarily designed to power the mp6-AIR pumps and allowing serial connection for enhanced pressure generation.

Flow rate and pressure generation of mp6-AIR, mp6-pp and mp6 depend on amplitude as well as on frequency. With gases the peak value for flow rate and pressure generation is at a higher frequency (first local maximum is at ~1250 Hz) than for liquids. Our standard electronic products have a maximal frequency of 226 Hz and 300 Hz, mp6-OEM and mp-x respectively, as they were designed to cover the pumping optimum for liquids at around 100 Hz. With the new mp6-QuadEVA the maximal frequency is now raised to 800 Hz, which allows much higher flow rates and pressure levels for gases.

One other difference to the standard electronics is that the mp6-QuadEVA generates a sine driving signal only. It is not as noisy and loud as for instance a rectangular signal. The evaluation board allows controlling


Technical data of mp6-QuadEVA	
Dimensions	80 x 60 x 16 cm 3.15 x 2.36 x 0.63 in.
Pumping media	liquids, gases
Max. volume flow (mp6-AIR, gas: air)	each single mp6-AIR: ~42 ml/min @ 260 V; 800 Hz
Adjustable parameters	Amplitude, frequency
Amplitude range	0 – 260 Vpp
Frequency range	50 – 800 Hz
Signal form	sine
Power supply	7.5 V, 1000 mA
Current consumption	avg. 220 mA, peak 280 mA ¹
¹ four connected mp6-AIR @ 260 V and 800 Hz	

amplitude and frequency either locally on board or via USB-connection. The local control is realized with rotary buttons with preset values for its positions. Each pump stage has its own on/off switch and status LED. The remote control via USB-connection can be achieved with any terminal software that is capable of sending commands to serial COM-Ports (e.g. HyperTerminal or PuTTY (www.putty.org)) but also with the software that comes together with the board. The board requires its own power supply regardless of local or remote control.

The evaluation board mp6-QuadEVA allows customers to experiment with multiple pumps for the achievability of their desired flow and pressure goals. Although gas applications are the main focus here, liquids can be pumped nevertheless.

For pumping air, the flow rate and pressure generation of a single mp6-AIR is displayed below. Additionally the forward pressure of serial connected mp6-AIR pumps is presented with the bottom right diagram. The cause for the pressure decrease of these lines is described further below.

Bartels Mikrotechnik offers the evaluation board as a set consisting of the mp6-QuadEVA board, a power supply, a mini-USB cable, four micropumps mp6-AIR together with software and USB driver.


Due to the character of an evaluation board concept the design and layout of the board is not optimized for any specific application or any serial device solution. The circuit can be adapted for various applications. On request we offer such an adaptation for your customized solution. Additionally it is possible to order the schematics, part list and coding for some license fee.

Why more frequency and why more pumps when pumping gas?


Especially for gas applications pressure levels are a key factor. The micropump mp6 was designed handling liquids, so its flow rate and pressure generation is best at a peak frequency of approximately 100 Hz (for water). With other frequencies the system "micropump mp6 with liquid inside" becomes less efficient and flow rate and pressure generation decreases.


Our micropumps mp6 and mp6-pp are measured with water at 100 Hz; the mp6-AIR is measured with air at 300 Hz. The reason here is that 300 Hz is the maximal frequency of our lab electronic mp-x.

More frequency and hence more performance of mp6-AIR is only possible with specialized electronic like the mp6-QuadEVA.



The diagrams below show flow rate and pressure generation – forward and suction pressure – and their dependency on frequency, different driving signals and amplitudes. Please note that the data is based on a waveform generator with an amplifier electronic; each data point represents the mean of five pumps.

The mp6-QuadEVA would achieve similar values, but due to an amplitude drop at higher frequencies and/or with higher load, i.e. connected pumps, the reached values are slightly lower. This amplitude behavior of the mp6-QuadEVA can be seen below (left hand picture) which is the reason for the pressure decrease mentioned above (and here with right hand picture).


So with higher frequencies more performance is possible but it may not be enough for some applications. Additionally, higher frequencies may be too noisy and loud or the increased power demand for driving higher frequencies is not available by the power source.

Fortunately, there exists an additional solution that can be used with or without higher frequency.

As mentioned above connecting pumps in serial succession increases the pressure level, this applies to liquids and gases. If all pumps are driven with the same signal and connected in a series the pressure generation adds up. Combining pumps in that way enhances forward and suction pressure.

This method cannot be continued forever as the actuators inside the pump have to be able to move under pressure load. When certain levels are passed the actuator stroke is decreasing hence flow rate and pressure generation too. This decrease begins roughly at 600 mbar for suction pressure and around 1000 mbar for forward pressure, so that further pumps will not increase the pressure linear anymore. The difference between suction and forward pressure is caused by the mechanical tensions of the actuator when bulging inwards or outwards.

Realizing such a serial connection of pumps with the mp-x (see diagram) requires one electronic for one pump. Using the mp6-QuadEVA will simplify the electronic side.

Pressure gain with pumps in serial connection, here with one mp-x for one mp6 pumping air.

Bartels Mikrotechnik can assist with additional engineering services bringing this evaluation board technology into your specific application. We can either offer to review your existing design for optimization or by a new development.

We are looking forward to discuss your individual needs in order to provide solutions tailored to your requirements. Contact us today!

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