

LABPLAS

Made in Canada

Lab Consumable Line

Sterile Sampling Bag

The Solution to Save Lab Samples



Food



Agriculture



Environment



Veterinary



Research



Pharmaceutical

REVODIX
레보딕스(주)

PRODUCT BRANDS

High-precision sampling for your industry

- Lab Plas는 다양한 실험 샘플을 보관 및 운송할 수 있는 가장 경제적이고 효율적인 방법을 제공합니다. 환경, 제약 연구, 품질 보증, 농식품 및 수의학 분야에서 사용 가능합니다.
- Lab Plas의 멸균 샘플링백은 누구나 간편하게 사용할 수 있습니다.
- Lab Plas의 멸균 샘플링백은 유리병에 비해 보관이 간편하며 샘플을 손쉽게 운송할 수 있습니다.

Pharmaceutical samples

- Lab Plas의 멸균 샘플링백을 이용하여 클린룸 장비를 멸균 상태로 이동 가능
- 클린룸에서 사용되는 원자재 보관 가능

Research samples

- 박테리아를 식별 가능
- 샘플을 혼합하는 과정에서 오염 방지
- 표본을 채취하거나 샘플 보관 가능

Food samples

- 교차오염을 방지하여 식품 분석 결과 신뢰도를 높임
- 원료를 효율적으로 보관하여 품질 관리 수준을 높임



Agricultural samples

- 토양 분석 샘플을 효과적으로 보관 가능
- U.S. Environmental Protection Agency (EPA) 인증을 취득

Veterinary samples

- 유전자 차이를 식별하기 위한 DNA 검사에 활용 가능
- 가축의 질병과 발생 가능한 질병의 연구에 활용 가능

TWIRL'EM BAGS

TWIRL'EM™ Sterile Sampling Bags

Labplas TWIRL'EM Sterile Sampling Bags은 액체, 고체 등 다양한 샘플을 보관할 수 있습니다. 외부 오염물질로부터 샘플을 보호하며 소프트한 외관을 가지고 있기 때문에 다양한 샘플을 보관할 수 있습니다.

- RNase, DNase, pyrogens, BPA Free



Solid

Semi-solid

Liquid


STERILITY:

- Extrusion process at high temperatures between 220°C and 240°C guarantees internal sterility;
- Product sterility verified through testing of random samples;
- In addition, sterilization by the extrusion process has been scientifically validated;
- R-NASE, D-NASE, Pyrogen and BPA free.

PRODUCT ADVANTAGES:

- Maximum bag mouth opening to facilitate insertion of sample;
- Lot number traceability for every product;
- Product meets USDA and ACIA industry standards and are made with FDA-compliant material.

Tear-off barrier
Ensure & maintain
sterility until usage of the bag.




Safety tabs. Ask for our safety tabs to avoid any possible puncture of the sample bag.



Choose from different combinations of wires to add resistance to your closing system.

- 2 round wires;
- 1 round wire & 1 flat wire;
- 2 flat wires.



Pull tabs



Write-on area available

STEP BY STEP INSTRUCTION

- 1 점선을 따라 샘플백의 상단을 찢어주세요
- 2 당김탭을 이용해 샘플백을 열어줍니다.
- 3 샘플을 샘플백에 주입합니다.
- 4 와이어의 양끝을 잡아 당겨 입구를 닫고 그대로 돌려줍니다.
- 5 와이어의 양쪽끝을 접어 단단히 고정하면 완벽하게 밀봉됩니다.

TWIRL'EM Sterile Sampling Bag

TWIRL'EM Sterile Sampling Bag



Safety Tabs – Clear – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
EPL-3050	500ea x 2pk	60	76 x 127	63
EPL-3070	500ea x 2pk	130	76 x 178	63
EPL-3570	500ea x 2pk	210	89 x 178	76
EPL-4575	500ea x 2pk	355	114 x 190	63
EPL-4590	500ea x 2pk	450	114 x 229	63
EPL-4512	500ea x 2pk	800	114 x 305	63
EPL-4515	500ea x 2pk	1200	114 x 382	63
EPL-5590	500ea x 2pk	650	140 x 229	76
EPL45590	500ea x 2pk	650	140 x 229	101
EPL-5515	500ea x 2pk	1500	140 x 382	76
EPL-7012	250ea x 4pk	1650	178 x 305	76
EPL47012	250ea x 4pk	1650	178 x 305	101
EPL-7015	250ea x 4pk	2500	178 x 305	76

Safety Tabs – Printed – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
EPR-3050	500ea x 2pk	60	76 x 127	63
EPR-3070	500ea x 2pk	130	76 x 178	63
EPR-3570	500ea x 2pk	210	89 x 178	76
EPR-4575	500ea x 2pk	355	114 x 190	63
EPR-4590	500ea x 2pk	450	114 x 229	63
EPR-4512	500ea x 2pk	800	114 x 305	63
EPR-4515	500ea x 2pk	1200	114 x 382	63
EPR-5590	500ea x 2pk	650	140 x 229	76
EPR45590	500ea x 2pk	650	140 x 229	101
EPR-5515	500ea x 2pk	1500	140 x 382	76
EPR-7012	250ea x 4pk	1650	178 x 305	76
EPR47012	250ea x 4pk	1650	178 x 305	101
EPR-7015	250ea x 4pk	2500	178 x 305	76

Safety Tabs – Clear – Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	µm
EFL-3070	500ea x 2pk	130	76 x 178	63
EFL-4590	500ea x 2pk	450	114 x 229	63
EFL-5590	500ea x 2pk	650	140 x 229	76
EFL-5512	500ea x 2pk	1080	140 x 305	76
EFL-7012	250ea x 4pk	1650	178 x 305	76
EFL47012	250ea x 4pk	1650	178 x 305	101
EFL-1012	250ea x 4pk	2250	254 x 305	101
EFL-1015	250ea x 4pk	4000	254 x 382	101

TWIRL'EM Sterile Sampling Bag

TWIRL'EM Sterile Sampling Bag



Safety Tabs – Printed – Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	µm
EFR-3070	500ea x 2pk	130	76 x 178	63
EFR-4590	500ea x 2pk	450	114 x 229	63
EFR-5590	500ea x 2pk	650	140 x 229	76
EFR-5512	500ea x 2pk	1080	140 x 305	76
EFR-7012	250ea x 4pk	1650	178 x 305	76
EFR47012	250ea x 4pk	1650	178 x 305	101
EFR-1012	250ea x 4pk	2250	254 x 305	101
EFR-1015	250ea x 4pk	4000	254 x 382	101

Regular Tabs – Clear – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
BPL-5590	500ea x 2pk	650	140 x 229	76
BPL-5512	500ea x 2pk	1080	140 x 305	76
BPL-5515	500ea x 2pk	1500	140 x 382	76
BPL-7012	250ea x 4pk	1650	178 x 305	76

Regular Tabs – Printed – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
BPR-5590	500ea x 2pk	650	140 x 229	76
BPR-5512	500ea x 2pk	1080	140 x 305	76
BPR-5515	500ea x 2pk	1500	140 x 382	76
BPR-7012	250ea x 4pk	1650	178 x 305	76

Regular Tabs – Clear – Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	µm
BFL-5590	500ea x 2pk	650	140 x 229	76
BFL45590	500ea x 2pk	650	140 x 229	101
BFL-5512	500ea x 2pk	1080	140 x 305	76
BFL-5515	500ea x 2pk	1500	140 x 382	76
BFL-7012	250ea x 4pk	1650	178 x 305	76

Regular Tabs – Printed – Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	µm
BFR-5590	500ea x 2pk	650	140 x 229	76
BFR-5512	500ea x 2pk	1080	140 x 305	76
BFR-5515	500ea x 2pk	1500	140 x 382	76
BFR-7012	250ea x 4pk	1650	178 x 305	76



TWIRL'EM Large Format



TWIRL'EM Large Format



Safety Tabs – Clear – Closure with 2 flat wires

P/N	Package	ml	mm x mm	µm
EDL41212	250ea x 1pk	3,000	305 x 305	101
EDL41214	250ea x 1pk	4,250	305 x 356	101
EDL41216	250ea x 1pk	5,500	305 x 406	101
EDL41218	250ea x 1pk	6,750	305 x 457	101
EDL41220	250ea x 1pk	8,000	305 x 508	101
EDL41224	250ea x 1pk	10,500	305 x 610	101
EDL41228	250ea x 1pk	13,000	305 x 711	101
EDL41520	250ea x 1pk	12,000	382 x 508	101
EDL41524	250ea x 1pk	15,000	382 x 610	101

Safety Tabs – Printed – Closure with 2 flat wires

P/N	Package	ml	mm x mm	µm
EDR41212	250ea x 1pk	3,000	305 x 305	101
EDR41214	250ea x 1pk	4,250	305 x 356	101
EDR41216	250ea x 1pk	5,500	305 x 406	101
EDR41218	250ea x 1pk	6,750	305 x 457	101
EDR41220	250ea x 1pk	8,000	305 x 508	101
EDR41224	250ea x 1pk	10,500	305 x 610	101
EDR41228	250ea x 1pk	13,000	305 x 711	101
EDR41520	250ea x 1pk	12,000	382 x 508	101
EDR41524	250ea x 1pk	15,000	382 x 610	101

Tear-off opening – Clear

P/N	Package	ml	mm x mm	µm
SCL21520	250ea x 2pk	12,000	381 x 508	51
SCL22424	250ea x 2pk	48,000	610 x 610	51
SCL33030	250ea x 1pk	55,000	760 x 760	76
SCL33037	250ea x 1pk	61,000	760 x 940	76
SCL43048	125ea x 1pk	81,000	760 x 1219	101



TWIRL'EM ECOLO

TWIRL'EM ECOLO



● ● ● Ecolo 제품은 100% 생분해성 LDPE으로 만들어진 환경 친화적인 제품입니다.
 Lab Plas 제품 중 녹색 스트랩이 둘러진 제품은 환경 친화적인 제품군 입니다.

Safety Tabs - Clear - Closure with 2 round wires

P/N	Package	ml	mm x mm	μm
EPL-3070E	500ea x 2pk	130	76 x 178	63
EPL-4590E	500ea x 2pk	450	114 x 229	63
EPL-5590E	500ea x 2pk	650	140 x 229	76
EPL45590E	500ea x 2pk	650	140 x 229	101
EPL-5515E	500ea x 2pk	1500	140 x 382	76
EPL-7012E	250ea x 4pk	1650	178 x 305	76
EPL47012E	250ea x 4pk	1650	178 x 305	101
EPL-7015E	250ea x 4pk	2500	178 x 382	76

Safety Tabs - Printed - Closure with 2 round wires

P/N	Package	ml	mm x mm	μm
EPR-3050E	500ea x 2pk	60	76 x 127	63
EPR-3070E	500ea x 2pk	130	76 x 178	63
EPR-4590E	500ea x 2pk	450	114 x 229	63
EPR-5590E	500ea x 2pk	650	140 x 229	76
EPR45590E	500ea x 2pk	650	140 x 229	101
EPR-7012E	250ea x 4pk	1650	178 x 305	76
EPR47012E	250ea x 4pk	1650	178 x 305	101

Safety Tabs - Clear - Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	μm
EFL-3070E	500ea x 2pk	130	76 x 178	63
EFL-4590E	500ea x 2pk	450	114 x 229	63
EFL-5590E	500ea x 2pk	650	140 x 229	76
EFL-5512E	500ea x 2pk	1080	140 x 305	76
EFL-7012E	250ea x 4pk	1650	178 x 305	76
EFL47012E	250ea x 4pk	1650	178 x 305	101
EFL-1012E	250ea x 4pk	2250	254 x 305	101
EFL-1015E	250ea x 4pk	4000	254 x 382	101

Safety Tabs - Printed - Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	μm
EFR-3070E	500ea x 2pk	130	76 x 178	63
EFR-4590E	500ea x 2pk	450	114 x 229	63
EFR-5590E	500ea x 2pk	650	140 x 229	76
EFR-7012E	250ea x 4pk	1650	178 x 305	76
EFR47012E	250ea x 4pk	1650	178 x 305	101
EFR-1012E	250ea x 4pk	2250	254 x 305	101
EFR-1015E	250ea x 4pk	4000	254 x 382	101



TWIRL' BLUE TWIRL'EM OPAQUE

TWIRL' BLUE

Polyethylene blue sterile bag

P/N	Package	ml	mm x mm	μm
EPR-4590-B	500ea x 2pk	450	114 x 229	63
EPR-5590-B	500ea x 2pk	650	140 x 229	76
EPR-7012-B	250ea x 4pk	1650	178 x 305	76



TWIRL'EM OPAQUE

Safety Tabs – Opaque – Closure with 2 round wires

P/N	Package	ml	mm x mm	μm
EPN-4590	500ea x 2pk	450	114 x 229	63
EPN-4512	500ea x 2pk	900	114 x 305	63
EPN-4515	500ea x 2pk	1200	140 x 382	63



TWIRL'EM CLEANROOM WATER SAMPLING KIT

TWIRL'EM Cleanroom

  **DOUBLE PACKAGING FOR INTERNAL AND EXTERNAL STERILITY!**

Safety Tabs – Clear – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
EPL-4512S	1000ea x 1pk	800	114 x 305	63
EPL-5590S	1000ea x 1pk	650	140 x 229	76
EPL-5515S	1000ea x 1pk	1500	140 x 382	76
EPL-7012S	1000ea x 1pk	1650	178 x 305	76

Safety Tabs – Printed – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
EPR-4512S	1000ea x 1pk	800	114 x 305	63
EPR-5590S	1000ea x 1pk	650	140 x 229	76
EPR-5515S	1000ea x 1pk	1500	140 x 382	76
EPR-7012S	1000ea x 1pk	1650	178 x 305	76

Safety Tabs – Clear – Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	µm
EFL-5590S	1000ea x 1pk	650	140 x 229	76
EFL45590S	1000ea x 1pk	650	140 x 229	101
EFL-5512S	1000ea x 1pk	1080	140 x 305	76
EFL-7012S	1000ea x 1pk	1650	178 x 305	76
EFL47012S	1000ea x 1pk	1650	178 x 305	101
EFL-1012S	1000ea x 1pk	2250	254 x 305	101
EFL-1015S	1000ea x 1pk	4000	254 x 382	101

Safety Tabs – Printed – Closure with 1 round wire & 1 flat wire

P/N	Package	ml	mm x mm	µm
EFR-5590S	1000ea x 1pk	650	140 x 229	76
EFR45590S	1000ea x 1pk	650	140 x 229	101
EFR-5512S	1000ea x 1pk	1080	140 x 305	76
EFR-7012S	1000ea x 1pk	1650	178 x 305	76
EFR47012S	1000ea x 1pk	1650	178 x 305	101
EFR-1012S	1000ea x 1pk	2250	254 x 305	101
EFR-1015S	1000ea x 1pk	4000	254 x 382	101

Water Sampling Kit

Safety Tabs – Stand-up bag – Printed – Closure with 2 round wires

P/N	Package	ml	mm x mm	µm
KWS-21100	100ea x 1pk	100	76 x 178	63
KWS-22200	100ea x 1pk	300	114 x 229	63



FILTRA-BAG

FILTRA-BAG

FILTRA-BAG with Open Top – Clear

P/N	Package	Description	ml	mm x mm	μm
SCT07012A	100ea x 4pk	For 400 series blenders	1650	178 x 305	101

FILTRA-BAG with Tear-Off Protection Strip – Clear

P/N	Package	Description	ml	mm x mm	μm
SCT-6090A	100ea x 4pk	For 400 series blenders	710	152 x 229	101
SCT-7012A	100ea x 4pk	For 400 series blenders	1650	178 x 305	101
SCT-1015A	100ea x 4pk	For 3500 series blenders	2720	254 x 381	101
SCT-1515A	100ea x 2pk	For 3500 series blenders	4000	381 x 381	101
SCT-1520A	100ea x 2pk	For 3500 series blenders	12 L	381 x 508	101

FILTRA-BAG – Safety Tabs – Clear – Closure with 1 round wire & 1 flat wire

P/N	Package	Description	ml	mm x mm	μm
EFT-6090A	100ea x 4pk	For 400 series blenders	710	152 x 229	101
EFT-7012A	100ea x 4pk	For 400 series blenders	1650	178 x 305	101
EFT-1015A	100ea x 4pk	For 3500 series blenders	2720	254 x 381	101
EFT-1515A	100ea x 2pk	For 3500 series blenders	4000	381 x 381	101
EFT-1520A	100ea x 2pk	For 3500 series blenders	12 L	381 x 508	101

STEP BY STEP INSTRUCTION

1 고품의 샘플을 필터백의 프린트가 있는 면으로 넣습니다.



2 균질기에 샘플백을 넣고 균질화 작업을 진행합니다.



3 내부필터를 거친 샘플을 피펫으로 추출합니다.





SANI-SPONGE

SANI-SPONGE



SANI-SPONGE Kit with Dry Sponge

P/N	Package	Description	ml	mm x mm	μm
KSS-61100	100ea x 4pk	Dehydrated sponge	450	114 x 229	63
KSS-61105	100ea x 4pk	Dehydrated sponge with gloves	450	114 x 229	63

SANI-SPONGE Kit with Wet Sponge


P/N	Package	Description	ml	mm x mm	μm
KSS-67110-DE	20ea x 5pk	DE-Neutralizing buffer	600	160 x 229	101
KSS-67110-LT	20ea x 5pk	Lethen broth	600	160 x 229	101
KSS-67110-NE	20ea x 5pk	Neutralizing buffer	600	160 x 229	101
KSS-67110-BPW	20ea x 5pk	Buffered peptone water	600	160 x 229	101
KSS-67115-DE	20ea x 5pk	DE-Neutralizing buffer with gloves	600	160 x 229	101
KSS-67115-LT	20ea x 5pk	Lethen broth with gloves	600	160 x 229	101
KSS-67115-NE	20ea x 5pk	Neutralizing buffer with gloves	600	160 x 229	101
KSS-67115-BPW	20ea x 5pk	Buffered peptone water with gloves	600	160 x 229	101

SANI-SPONGE sampling kits are produced in accordance with HACCP, USDA, ACIA and CFIA requirements




STEP BY STEP INSTRUCTION


1 SANI-SPONGE의 겉면 포장을 조심히 뜯어 제품을 꺼냅니다. (글러브 착용을 권장합니다.)




2 스폰지의 한쪽면으로 표면 샘플링을 진행합니다. (30*30cm 의 구역을 권장합니다.)




3 스폰지의 반대면으로 표면 샘플링을 반복합니다.



4 샘플링이 끝난 스폰지를 다시 샘플링 백에 넣고 입구를 봉합니다.



5 제품 겉면에 샘플에 대한 정보를 기입합니다.



SANI-STICK

SANI-STICK



SANI-STICK Kit with Wet Sponge

P/N	Package	Description	ml	mm x mm	μm
KSS-67310-DE	10ea x 10pk	Sponge on handle with DE-Neutralizing buffer	600	160 x 229	101
KSS-67310-LT	10ea x 10pk	Sponge on handle with Lethen broth	600	160 x 229	101
KSS-67310-NE	10ea x 10pk	Sponge on handle with Neutralizing buffer	600	160 x 229	101
KSS-67310-BPW	10ea x 10pk	Sponge on handle with Buffered peptone water	600	160 x 229	101
KSS-67315-DE	10ea x 10pk	Sponge on handle with DE-Neutralizing buffer with gloves	600	160 x 229	101
KSS-67315-LT	10ea x 10pk	Sponge on handle with Lethen broth with gloves	600	160 x 229	101
KSS-67315-NE	10ea x 10pk	Sponge on handle with Neutralizing buffer with gloves	600	160 x 229	101
KSS-67315-BPW	10ea x 10pk	Sponge on handle with Buffered peptone water with gloves	600	160 x 229	101

Simply squeeze and release the sponge



1.5" x 3" biocide-free cellulose sponge gamma irradiated. Pre-moistened with a variety of broths or dry



Polypropylene rigid handle for sampling hard-to-reach areas



No plastic residue will remain on sponge after sampling.
Friendly to blending procedures.

STEP BY STEP INSTRUCTION

1. 제품의 포장을 조심히 제거하고 SANI-STICK의 손잡이를 잡고 꺼냅니다.
2. 젖은 스폰지를 이용해 표면을 닦아냅니다.
3. 샘플백에 스폰지를 넣고 손잡이의 끝을 잡아당기면 스폰지가 분리됩니다.
4. 멸균백을 3~4회 정도 감아 밀봉합니다.
5. 샘플의 정보를 기입합니다.

SECURE-T SECURE-STRIP



SECURE-T™

SECURE-T with Tear-Off Protection Strip – Clear

P/N	Package	Description	ml	mm x mm	µm
SCL-4060	250ea x 4pk	For 80 series blenders	205	102 x 152	76
SCL-5512	500ea x 2pk	For 400 series blenders	1080	140 x 305	76
SCL-7012	250ea x 4pk	For 400 series blenders	1650	178 x 305	76
SCL47012	250ea x 4pk	For 400 series blenders	1650	178 x 305	101
SCL41214	250ea x 2pk	For 3500 series blenders	4.5 L	305 x 356	101
SCL41218	250ea x 1pk	For 3500 series blenders	6.75 L	305 x 457	101
SCL-1520	250ea x 2pk	For 3500 series blenders	12 L	381 x 508	76
SCL41520	250ea x 2pk	For 3500 series blenders	12 L	381 x 508	101

SECURE-T with Tear-Off Protection Strip – Printed

P/N	Package	Description	ml	mm x mm	µm
SCR-7012	250ea x 4pk	For 400 series blenders	1650	178 x 305	76

BLENDER BAG with Open Top – Clear

P/N	Package	Description	ml	mm x mm	µm
SCL07012	250ea x 4pk	For 80 series blenders	1650	178 x 305	76
SCL01520	250ea x 2pk	For 400 series blenders	12 L	381 x 508	76

BLENDER BAG with Open Top – Printed

P/N	Package	Description	ml	mm x mm	µm
SCR07012	250ea x 4pk	For 80 series blenders	1650	178 x 305	76

SECURE-STRIP

P/N	Package	Description	ml	mm x mm	µm
SCLR7012-B	Box dispenser	For 400 series blenders	1650	178 x 305	76
SCLR7012	4 rolls of 250 units	For 400 series blenders	1650	178 x 305	76



Looking for a permanent solution? Ask us about our stainless steel dispenser. Another way to make blender bags easy to use in the workplace!



SAMPLING ACCESSORIES TO MAKE LAB WORK EASIER!
Labplas offers a comprehensive range of accessories designed to simplify sample processing and increase efficiency in the lab.

CUSTOM BAGS

Customization and Private Branding



자유로운 프린팅: 로고, 이미지 추가 가능

Create brand awareness by printing your logo on the bags!
Let us know your specifics and we will create a unique bag for your needs

- Two-color limitation.

원하는 사이즈로 제작 가능

- Thickness of 2.5, 3.0, 3.5 or 4 mil.in or more;
- Minimum width of 3 inches, no maximum;
- Minimum length of 5 inches, no maximum.



자유로운 패키징 가능: 패키징 당 수량 선택 가능, 클린룸 전용 패키징 가능

- Units per box ;
- Clean room special packaging ;
- Units of sponges or special inserts per bag.

Minimum quantity to be determined upon quote of product.

KIT 제작 가능: 멸균 제품 포함하여 제작 가능 (Spoon, swabs, clippers 등)

Labplas offers the option of customized insertions such as:

- Wet sponges (your unique formula of media);
- Sterile objects such as spoons, swabs, clippers, etc.;
- Pre-rolled bags for diagnostic kits;
- Cotton tissue with several choices of media for analysis;
- Boot swabs with several choices of media for analysis;
- Any other object you would like to be inserted in your bags.

ACCESSORIES

Rack'N'Lab

Product code: SA5-3070, SA5-7012



Labrack

Product code: SA4-7012



Closure Clips

Product code: SA2-4060 (80 SERIES)
SA2-7012 (400 SERIES)
SA2-1520 (3500 SERIES)



Work Station

Product code: SA1-7012



Gloves

Product code: SA7-GLOVES



Pen

Product code: SA6-DETECPEN





TWIRL'EM[®]

TWIRL'EM[®]
ECOLO

TWIRL'BLUE

SECURE-T[®]

FILTRA-BAG[®]

LAB/RACK

Sani-STICK[®]

SANI SPONGE[®]

Rack Lab

Contact us for all your sampling needs.

