



Dosage of Liquids: Series T

1. Zero Dead Volume Syringes

Plunger: Stainless steel, plunger runs inside the needle
Glass Barrel: Borosilicate 3.3
Precision: $\pm 1\%$ of the volume
Main Use: - sample dosage for gas chromatography
 - dosage of small amounts of liquids
 - for Thin Layer Chromatography (TLC)

Volume	Tip Style Length O.D.	a 70 mm	b 70 mm	d 70 mm
1 μ l	0.47 mm	2200101	2200102	2200103
5 μ l	0.55 mm	2200301	2200302	2200303

Volume	Tip Style Length O.D.	a 80 mm	b 80 mm	d 80 mm
2 μ l	0.72 mm	2200201	2200202	2200203
10 μ l	0.70 mm	2200401	2200402	2200403
25 μ l	1.10 mm	2200501	2200502	2200503

Dosage of Liquids: Series T

2. Replacement Parts

Plunger: Stainless steel, plunger runs inside the needle
Glass Barrel: Borosilicate 3.3
Precision: $\pm 1\%$ of the volume
Main Use: - sample dosage for gas chromatography
 - dosage of small amounts of liquids
 - for Thin Layer Chromatography (TLC)

Replacement Kit

Includes 1 needle and 1 calibrated plunger

Volume	Tip Style Length O.D.	a 70 mm Needle/ Plunger	a 80 mm Needle/ Plunger	b 70 mm Needle/ Plunger	b 80 mm Needle/ Plunger
1 μ l	0.47 mm	2231101	-	2231102	-
2 μ l	0.72 mm	-	2231201	-	2231202
5 μ l	0.55 mm	2231301	-	2231302	-
10 μ l	0.70 mm	-	2231401	-	2231402
25 μ l	1.10 mm	-	2231501	-	2231502

Sealing Kit

Includes 5 seals

Volume	O.D.	
1 μ l	0.47 mm	2230001
2 μ l	0.72 mm	2230002
5 μ l	0.55 mm	2230003
10 μ l	0.70 mm	2230004
25 μ l	1.10 mm	2230005